Application and safety instructions for the MEISINGER Periimplantitis Kit

CONTENT

Language	Page
English (EN)	2-4
French (FR)	5-7
German (DE)	8-10
Italian (IT)	11-13
Spanish (ES)	14-16
Turkish (TR)	17-19

Application and safety instructions for MEISINGER

Periimplantitis Kit

Art.-No.: 2575

The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the ideal peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counterclockwise so the implant surface is smoothed.

INDICATIONS FOR USE

The Bone Management Periimplantitis Kit is indicated for the intraoral processing and smoothing of exposed titanium implant surfaces for the surgical treatment of peri-implantitis.

UTILIZATION

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Preparation

The initial situation is an exposed titanium implant surface after bone resorption due to existing peri-implantitis. As preparation for the optimal peri-implantitis treatment, the prosthetics have to be removed from the implant. The implant is provided with a native cover screw, so that the interface is protected against contamination.

Processing of the titanium implant surface

Prior to the processing of the titanium surface, the contaminated tissue around the implant is removed with the aid of an appropriate manual instrument (e.g. Curette Art.-No. HI317). Then, the surface structure of the implant is removed with the fine tungsten carbide finishers before the titanium is smoothed with the ultra-fine finishers. During the processing **the clockwise rotating instruments are led around the exposed surface counterclockwise** and with low contact pressure, without damaging the implant shoulder. This kind of use allows for the best possible control of the used instruments. Size and shape of the instruments depend on the existing patient situation and the implant shape.

When using the ultra-fine finishing bur a very smooth surface structure is achieved. The use of elastic polishers can and should be avoided, also to avoid silicone residues.

Post-processing

Prior to the tension-free wound closure all metal residues have to be removed by sufficient rinse and/or the application of appropriate prophylaxis brushes (e.g. Art.-No. 115, 116).

Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. Dual Waste System).

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Adverse reaction

By smoothing the implant surface, contamination with metal shavings can not be completely avoided. This may possibly promote wound healing disorders.

NOTE: This procedure is a surgical therapy. Therefore, all alternative, non-surgical treatment options should be considered in advance and the best possible procedure, according to available data, should be selected. If the progression of the inflammation has already progressed too far or should not be stopped in the treatment, it should be assessed whether an explantation of the implant should be performed in order to avoid further bone resorption. During the whole processing time, attention has to be paid to sufficient cooling with sterile cooling liquid. Furthermore, the implant shoulder may not be damaged, so that the implant-abutment-interface remains intact also after the treatment. Excessive pressure during the use of the instruments has to be avoided. Also the use of the instruments with leverage effect is not allowed, due to an increased risk of breakage.

CAS: 7440-48-4
Contains hazardous
substances:
Cobalt,
CAS: 7440-48-4

This product contains more than 0.1 % by mass of cobalt and therefore must be labelled as CMR substance class 1B (carcinogenic, mutagenic, reprotoxic) according to current regulations. It has been demonstrated that there is no increased risk of cancer or adverse reproductive or genetic effects when used as directed.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Hager & Meisinger GmbH
Hansestr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

CONTENT

Tungsten carbide finishing burs

	1	1	1	1	1	1	1	1
Fig.	HM379	HM379	HM48L	HM48L	HM379U	HM379U	HM48LU	HM48LU
Shank ¹	317	317	317	317	317	317	317	317
Size ²	014	023	014	023	014	023	014	023
Length mm	3.1	4.2	8.0	8.0	3.1	4.2	8.0	8.0
	1.4	2.3	1.4	2.3	1.4	2.3	1.4	2.3
Opt. speed rpm	40.000	40.000	40.000	40.000	40.000	40.000	40.000	40.000
Max. speed rpm	190.000	190.000	120.000	120.000	190.000	190.000	120.000	120.000
CE 0044 / 								

¹ 317=FG XXL

² Largest working part diameter in 1/10 mm

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