Application and safety instructions for the MEISINGER

Drill-Stop-Control Drill-Stop-Control BLT Drill-Stop-Control myplant two

CONTENT

Language	Page
English (EN)	2-5
French (FR)	6-9
German (DE)	10-13
Italian (IT)	14-17
Spanish (ES)	18-21
Turkish (TR)	22-25

Application and safety instructions for MEISINGER

Drill-Stop-Control Drill-Stop-Control BLT Drill-Stop-Control myplant two

Art.-No.: BDS00, BDS02, BDSMP

The Bone Management Drill-Stop systems are used for simple and safe preparation of the implant site. The increasing diameters of the twist drills of the respective systems are perfectly matched. The stop sleeves are used for the defined depth stop during the pilot drilling as well as during the extension drilling. For use without stop sleeves, the instruments are equipped with laser markings.

INDICATIONS FOR USE

The Bone Management Systems Drill-Stop-Control, Drill-Stop-Control BLT and Drill-Stop-Control myplant two are designed for the surgical extension of the implant site in the upper and lower jaw in preparation for the insertion of a dental implant suitable for the respective system:

- **Drill-Stop-Control:** The diameters of the twist drills are matched to cylindrical implants of the Oktagon implant system and comparable implant systems
- **Drill-Stop-Control BLT:** The diameters of the twist drills are adapted to bone level tapered implants of the Oktagon implant system or comparable implant systems
- **Drill-Stop-Control myplant two:** The diameters of the twist drills are matched to the myplant two implant system

Please also observe the information provided by the implant manufacturer regarding indications and contraindications and the use of the selected implant system.


UTILIZATION

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Determination of the drilling depth / selection of the stop sleeve

The required drilling depth and selection of the appropriate stop sleeve is determined before the start of the operation according to the selected implant length. The specifications of all lengths of the stop sleeves contained in the system can be found in the article overview under “Content”.

⚠ **ATTENTION:** To avoid injuring anatomical structures, the apical overlength of the twist drills must be observed. You will find the relevant information in the article overview under “Content”.

Preparation of the implant site

For the subsequent preparation of the implant site, the selected drill sleeve must be placed on the drill used in each step. Make sure that the drill sleeve is securely seated. In order to avoid overheating of the jawbone, sufficient cooling with preferably physiological, sterile saline solution must be ensured during preparation. Furthermore, the drills are to be used intermittently and under low pressure.

Pilot drilling of the alveolar ridge

In order to prevent the twist drills from slipping when, for example, the alveolar ridge tapers to a point, the pilot drilling can be used beforehand.

Use of the extension drills

For the subsequent extension drilling, the twist drills are used in ascending diameter until the desired diameter of the implant site is reached. The first step is to drill a pilot hole with the smallest twist drill. This determines the alignment of the implant site in addition to the predefined depth. The final diameter as well as the fine preparation of the implant site must be created with the instruments of the planned implant system. Depending on the existing case situation and the desired implant diameter, the instruments to be used can be selected individually.

Implant insertion

Following the final preparation of the implant site with the system's own instruments, the implantation of the selected implant can be performed in accordance with the manufacturer's information.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical measures must be observed. In addition, the contraindications for the implantation of an endosseous dental implant should be considered. The respective implant manufacturer provides information about these in the instructions for use of the selected implant system. There are no specific contraindications for the use of the instruments contained in this kit.

 **IMPORTANT:** Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection / dehiscences).

 **ATTENTION:** The instruments are not intended for final implant site preparation. The diameter as well as the fine preparation of the implant site must be produced with the instruments of the planned implant system according to the manufacturer's instructions. Surgical and implantological procedure must always be performed according to the current state-of-the-art. This applies to all drilling operations: To reduce the risk of bone overheating and thus necrosis, the respective instrument must be used intermittently under low pressure and with constant cooling with sterile physiological saline solution.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disposal

When disposing of prosthetic components and instruments (at the end of their service life or after the specified shelf life has expired), care should be taken to ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national requirements (e.g. dual waste system).

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Hager & Meisinger GmbH
 Hansemannstr. 10
 41468 Neuss • Germany
 Tel. +49 2131 2012-0
 Fax: +49 2131 2012-222
 Email: info@meisinger.de
 Internet: www.meisinger.de

CONTENT





Drill-Stop-Control:

	Initial burr with stop	Twist Drills			
	1	1	1	1	1
Fig.	187RF	TDS15	TDS15	TDS15	TDS15
Shank ¹	204	204	204	204	204
Size ²	018	020	028	035	042
Length mm	19.0	19.0	19.0	19.0	19.0
	1.8	2.0	2.8	3.5	4.2
Opt. speed rpm	600-800	800-1.000	400-600	400-600	400-600
Max. speed rpm	1.000	1.500	800	800	800
CE 0044					

	Stop sleeves					
						
	1	1	1	1	1	1
Fig.	CL031	CL032	CL033	CL034	CL035	CL036
Length mm	4.0	5.0	6.0	7.0	8.0	9.0
≙ max. drilling depth mm	15.0	14.0	13.0	12.0	11.0	10.0
						

ATTENTION: To avoid injuring anatomical structures, the apical overlength of the twist drills of up to 0.8 mm must be observed.

Drill-Stop-Control BLT:

						
						
	1	1	1	1	1	1
Fig.	187RF	TDS0T	TDS0T	TDS0T	TDS0T	TDS0T
Shank ¹	204	204	204	204	204	204
Size ²	018	022	028	033	035	042
Length mm	19.0	19.0	19.0	19.0	19.0	19.0
	1.8	2.2	2.8	3.25	3.5	4.2
Opt. speed rpm	600-800	400-600	400-600	400-600	400-600	400-600
Max. speed rpm	1.000	800	800	800	800	800
CE 0044						
						
	1	1	1	1	1	1
Fig.	CL061	CL062	CL063	CL064	CL065	CL066
Length mm	5.0	7.0	9.0	11.0	13.0	15.0
$\triangleq$ max. drilling depth mm	14.0	12.0	10.0	8.0	6.0	4.0
CE						

ATTENTION: To avoid injuring anatomical structures, the apical overlength of the twist drills of up to 0.4 mm must be observed.

CONTENT

Drill-Stop-Control myplant two:

	Initial bur with stop		Twist Drills							
	1	1	1	1	1	1				
Fig.	187RF	TDS18	TDS18	TDS18	TDS18	TDS18				
Shank ¹	204	204	204	204	204	204				
Size ²	018	020	024	029	033					
Length mm	19.0	22.0	22.0	22.0	22.0					
	1.8	2.0	2.4	2.9	3.3					
Opt. speed rpm	600-800	800-1.000	600-800	400-600	400-600					
Max. speed rpm	1.000	1.500	1.000	800	800					
	CE 0044									
	Stop sleeves									
	1	1	1	1	1	1	1	1	1	1
Fig.	CL050	CL051	CL058	CL052	CL053	CL054	CL055	CL056	CL057	CL059
Length mm	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	14.0	15.0
max. drilling depth mm	16.0	15.0	14.0	13.0	12.0	11.0	10.0	9.0	8.0	7.0
	CE									

⚠ ATTENTION: To avoid injuring anatomical structures, the apical overlength of the twist drills of up to 0.6 mm must be observed.

¹ 204=RA ² Largest working part diameter in 1/10 mm

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