
# Application and safety instructions for the MEISINGER Trephine Ejection Kit

## CONTENT

| Language           | Page  |
|--------------------|-------|
| English (EN) ..... | 2-4   |
| French (FR) .....  | 5-7   |
| German (DE) .....  | 8-10  |
| Italian (IT) ..... | 11-13 |
| Spanish (ES) ..... | 14-16 |
| Turkish (TR) ..... | 17-19 |

## Application and safety instructions for the MEISINGER Trephine Ejection Kit

Art.-No.: BTE00

The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special: Two-part, internally cooled trephines allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trephines. Furthermore, the included pre-drilling trephines allow for a safe application of the trephines. Due to the two-part trephines, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

### INSTRUCTIONS FOR USE

The Bone Management system Trephine Ejection Kit is indicated for the intraoral extraction of bone cylinders for augmentation of bone deficits in the upper and lower jaw.

### INSTRUCTION

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.


#### Before Application

The instruments are supplied non-sterile, so they need to be cleaned, disinfected and sterilized prior to use, this applies in particular to the first-time use. For cleaning and disinfection the products are disassembled as far as possible. For that the shank of the trephine is attached to the hand wheel, so the trephine working part can be released and removed by a left-hand turn. The rubber ring remains on the shank during the first processing. After cleaning and disinfection have been completed, the trephine working part is attached again to the shank and fixed by a right-hand turn. The secure fit of the working part has to be checked. Subsequently, the instruments can be sterilized according to the instructions for processing. Since the trephines have an internal cooling, these are instruments with increased processing requirements. According to the RKI guideline, mechanical cleaning and disinfection in disinfection devices and sterilization with moist heat is recommended. For proper processing of the articles, please note the instructions for processing (cleaning, disinfection and sterilization) of medical devices produced by Hager & Meisinger GmbH.

#### Initial Drilling

According to the planned diameter of the bone cylinder the appropriate initial drill is selected. With this drill the bone at the donor site is punch-marked in order to avoid slippage of the trephine.


#### Extraction of the bone cylinder from the donor site

According to the planned diameter of the bone cylinder the appropriate trephine is selected. Before use, the fixed position of the trephine working part on the shank has to be checked.

**Caution: Depending on the construction, the tre-phines may not be used counterclockwise.**


At the punch-marked donor site the bone cylinder of the planned length can now be extracted with the aid of the trephine. For this purpose, the trephine is placed on the cortical bone and intermittently inserted into the bone. It is recommended to support the working hand to control the trephine optimally. The speed should not be too low to avoid tilting of the trephine in the depth of the bone. Please note the recommended and maximum speeds. The depth marks on the trephine working parts serve for length orientation.

#### Extraction of the bone cylinder from the trephine

To extract the bone cylinder from the trephine working part, the trephine is removed from the handpiece. To achieve a better handling, the ejection needle is attached to the hand wheel. Then, it can be led through the lumen of the trephines' internal cooling from the shank end in the direction of the trephine working part, so the bone cylinder is pushed out.

**Caution:** In the case of soft bone, the application of the ejection needle is contraindicated due to its small diameter. In this case, proceed according to step "Alternative Extraction of the bone cylinder from the trephine in case of soft bone".

#### Alternative Extraction of the bone cylinder from the trephine in case of soft bone

In case of soft bone, the extraction of the bone cylinder is performed with the aid of the ejection instruments. For that, the shank of the trephine is attached to the hand wheel, so that the trephine working part can be released and removed by a left-hand turn. Then, the trephine working part is inserted into the ejection sleeve. With the aid of the appropriate ejection instrument the bone cylinder is pushed out.

The 2.1 mm and 2.5 mm diameter trephines go with the 2.1 mm diameter ejection instrument. The 2.9 mm and 3.1 mm diameter trephines go with the 2.9 mm diameter ejection instrument.

#### After use

The instruments need to be cleaned, disinfected and sterilized after every application. For cleaning and disinfection the products are disassembled as far as possible. With the aid of a tweezer the rubber rings on the shank are now removed and disposed. Then, the instruments can be reprocessed. After completed cleaning and disinfection, preferably completely mechanical, a new rubber ring (Art.-No. 2151) is mounted on the shank and the trephine working part is again attached and fixed to the shank by a right-hand turn. Then, the instruments can be sterilized according to the instructions for processing. For a proper processing of the articles please note the instructions for processing (cleaning, disinfection and sterilization) of medical devices produced by Hager & Meisinger GmbH.

### Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

#### Absolute contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes
- Insufficient bone supply (quality / quantity) of the donor site as well as insufficient bone quality of the donor region

#### Relative contraindications

- Diseases that affect bone metabolism
- Drug and alcohol abuse
- Lack of cooperation of the patient
- Poor circulation
- Hard work or active sports
- Mental condition, which can lead to disregard of the medical order
- Highly atrophic jaw

**IMPORTANT:** Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection / dehiscences).

For all cutting and drilling procedures (transplant removal / other work), the following applies: In order to reduce the risk of overheating the bone and necrosis formation along with it, the respective instrument has to be used intermittently using little pressure and continuous cooling using a sterile physiological saline solution.

It should be avoided to trepan the lingual cortical bone as damage to the soft tissue and vessels of the mouth base can lead to serious complications. The extracted bone graft should be used as quickly as possible in the recipient site or stored in blood or isotonic saline solution. It is recommended to fill larger defects after bone extraction. In the donor site, a tension-free and complete wound closure must follow quickly.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national requirements (e.g. dual waste system).

### General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Manufacturer

REF

Item number

MD

Medical device

LOT

Batch number

UDI

Unique Device Identification

Do not re-use

Non-sterile

Rx only

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Caution

CE / CE 0044

Conformity mark

Do not use if the packaging is damaged

Do not use if the packaging is damaged

Manufacturing date

Consult electronic instructions for use

ifu.meisinger.de

**Hager & Meisinger GmbH**

Hansestr. 10  
41468 Neuss • Germany  
Tel: +49 2131 2012-0  
Fax: +49 2131 2012-222  
Email: [info@meisinger.de](mailto:info@meisinger.de)  
Internet: [www.meisinger.de](http://www.meisinger.de)

## CONTENT

|                    | Initial bur trephines |         |         |         | Shank   | Trephines |         |         |         |
|--------------------|-----------------------|---------|---------|---------|---------|-----------|---------|---------|---------|
|                    |                       |         |         |         |         |           |         |         |         |
|                    | I                     | I       | I       | I       | 4       | I         | I       | I       | I       |
| Fig.               | 230KH                 | 230KH   | 230KH   | 230KH   | 229KH*  | 229KH     | 229KH   | 229KH   | 229KH   |
| Shank <sup>1</sup> | 204                   | 204     | 204     | 204     | 204     | 900       | 900     | 900     | 900     |
| Size <sup>2</sup>  | 02I                   | 025     | 029     | 033     | 000     | 02I       | 025     | 029     | 033     |
| Length mm          | 2.0                   | 2.0     | 2.0     | 2.0     | 20.5    | 15.5      | 15.5    | 15.5    | 15.5    |
|                    | 2.1                   | 2.5     | 2.9     | 3.3     | -       | 2.1       | 2.5     | 2.9     | 3.3     |
|                    | 3.1                   | 3.5     | 3.9     | 4.3     | -       | 3.1       | 3.5     | 3.9     | 4.3     |
| Opt. speed rpm     | 300-500               | 300-500 | 300-500 | 300-500 | 300-500 | 300-500   | 300-500 | 300-500 | 300-500 |
| Max. speed rpm     | 800                   | 800     | 800     | 800     | 4.000   | 2.000     | 2.000   | 2.000   | 2.000   |
| CE 0044            |                       |         |         |         |         |           |         |         |         |

|                    | Ejection instruments |       |       | Ejection sleeve | Hand wheel |
|--------------------|----------------------|-------|-------|-----------------|------------|
|                    |                      |       |       |                 |            |
|                    | I                    | I     | I     | I               | I          |
| Fig.               | ASTKH                | ASTKH | ASNKH | ASHKH           | HRKH I     |
| Shank <sup>1</sup> | -                    | -     | -     | -               | -          |
| Size <sup>2</sup>  | 02I                  | 029   | 008   | -               | -          |
| Length mm          | 22.0                 | 22.0  | 37.0  | 28.0            | 15.5       |
|                    | -                    | -     | -     | 5.5             | -          |
|                    | 5.5                  | 5.5   | -     | -               | -          |
|                    | 2.1                  | 2.9   | 0.8   | -               | -          |
| CE                 |                      |       |       |                 |            |

<sup>1</sup> 204=RA <sup>2</sup> Largest working part diameter in 1/10 mm

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