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Application and safety instructions for the MEISINGER
Screw System TX
Screw System TX Professional

EN

Screw System for Bone Fixation

Art.-No.: BTX00, BTXPR


The Osteosynthesis Systems Screw System TX and Screw System TX Professional allow for the universal use for fixation of autologous bone cylinders, bone blocks and bone chips for the purpose of pre-implant augmentation. The bone screws consist of titanium alloy Grade 5. The Torx connection guarantees reliable accommodation of the screws in the driving tool. At the same time, high tensile forces can be transferred. For each 0.9 mm screw (TCT09) available, the kits contain a rescue screw with a larger diameter (TCT10) in case the transplant does not allow high tensile fixation. The small diameter and the low head height of the screws allow for a comfortable application for patient and user. The Screw System TX Professional contains additional screws in the lengths of 13.0 mm and 16.0 mm.

Indications for Use

The Bone Management® systems Screw System TX and Screw System TX Professional are indicated for the fixation of autologous bone cylinders, bone blocks and bone chips for the purpose of pre-implant vertical or horizontal augmentation.


Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.


Preparation of the donor and recipient site

To optimize the healing process, the recipient site is prepared first with small dotted drillings (Bleeding Points). The bone transplant is adjusted to the bone site, which has been appropriately prepared, so that it lies flat on the recipient bone. Good contact between bone transplant and recipient site as well as the mechanical stability and the absence of mechanical movement of the transplant are requirements for the success of the treatment.


Placing the drill hole through bone cylinder and recipient site

The apical end of the bone screws is non-self-drilling. For this reason, a pre-drilling is required for the application of bone screws. For this purpose, a 0.9 mm drill hole is placed into the bone cylinder using the 203RF 205 009 Twist Drill (or for rescue screws a 1.0 mm hole with TC203). The drill hole is made up to the opposing cortical bone of the recipient site.

Placing of the gliding hole

Then, with the aid of the Twist Drill 203RF 205 011 (or A2001 for rescue screws) solely the bone cylinder is provided with a gliding hole. The recipient site is not provided with a gliding hole in order to achieve in the next step a fixation of the screw in the bone according to the lag screw technique.

Fixation of the bone transplant

Using the screw driver, a screw of the appropriate length is selected and the bone transplant is attached to the recipient bone. The screw can be inserted into the bone without threading, as it has a self-tapping thread design. It should be gently tightened until the bone cylinder is fixed firmly. Doing this, the bone may not be overloaded. If a secure fixation of the transplant with only one screw is not possible, the fixation with at least 2 screws is recommended for rotational security.

Postoperative proceeding

Following clinically and radiographically confirmed, complete bone healing and functional restoration, the screws are either removed or left in situ, depending on the patient situation and treatment goal. If there is an infection or weakening of the bone tissue through the inserted bone screw, the removal of the screw is required. The removal of the screws is also necessary in case of screw loosening and always in the treatment of children, in order to avoid growth inhibition.

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Contraindications
Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications
• Active infection in or near the fixation site as well as local or systemic pathological processes (e.g. symptoms like fever, local inflammation, abscess)
• Insufficient bone quality and/or quantity
• Diseases or medication that impair bone metabolism

Relative Contraindications
• Dentoalveolar growth which has not come to an end (exception: Cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
• Known allergy and / or foreign body reaction to the implant material
• Systemic and metabolic disorders or medical treatments leading to progressive deterioration of bone (e.g. cortisone, immunosuppressive, Bisphosphonates, Diabetes)
• Drug and alcohol abuse
• Lack of compliance by the patient
• Poor circulation
• Heavy labor or active sports
• Mental conditions that may cause patient's failure to adhere to physician's order
• Highly atrophic jaw

IMPORTANT: Attention must be paid to protecting anatomical structures (safety distance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection / resorption). Care should be taken to avoid contact between the transplant and adjacent teeth to prevent the development of infections and dehiscence. Cases where the mechanical stability and the absence of mechanical movement of the transplant can not be achieved may affect successful bone healing.

SEPARATE INSTRUCTIONS FOR SCREWS: Bone screws are compatible with the Screw System TX and Screw System TX Professional instruments. The screws may only be used with these system's own instruments. When selecting the screw length and the screw placement, the evaluation of the available bone volume by means of X-ray diagnostics must be taken into account. The screw head is to be inserted as flat as possible in the bone block to prevent the soft tissue prophylactically from perforation or dehiscences. Excessive physical activity as well as trauma affecting the augmentation site can lead to premature screw failure due to loosening or breaking. The screws may only be used in accordance with the application and safety instructions. Any kind of modification in the procedure or on the product may impair the success of the treatment or damage the product. The screws are a disposable product and may therefore not be reprocessed and reused after use. Even with already used, but optically undamaged screws, it may lead to microscopic damage due to the previous mechanical stress and thus to fatigue and product failure. The screws have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the screws in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

General instructions
Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Hager & Meisinger GmbH
Hansemannstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

Meisinger USA, L.L.C.
10150 E. Easter Avenue
Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400
Fax: +1 (303) 268-5407
Email: info@meisingerusa.com
Internet: www.meisingerusa.com

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Inhalt / Content / Contenido / Contenu / Contenuto

Table with 6 columns for different drill types: Spiralbohrer für TCT09, Spiralbohrer für TCT10, Hartmetall Kugelfräser, and Schraubendreher. Includes technical drawings and a table with columns for Fig., Shank, Size, Length mm, Δ, Opt. speed rpm, and Max. speed rpm. CE 0044 logo at the bottom.

Zusätzlicher Inhalt / Additional Content / Contenido adicional / Contenu complémentaire / Contenuto aggiuntivo

Table for Screw System TX Professional showing additional content for screws. Columns include Fig., Shank, Size, Length mm, Δ, and Opt. speed rpm. CE 0044 logo at the bottom.