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With the aid of a fitting pilot bur (A2001 for the Extraction screws BE001230 060 and BE001230 070, B2001 for the Extraction screws BE001230 080 and BE001230 085) an axial drill hole has to be drilled into the center of the dental tooth root or tooth root fragment. This should have a depth of about 7 mm in hard tissue. A deeper drill hole is not necessary. The use is carried out with sufficient cooling, preferably with a physiological, sterile saline solution. In order to remove the drilling chips easily, the pilot bur has to be used intermittently, i.e. with up and down movements.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Preparation of the tooth root

The desmodontal fibers (Sharpey's fibers) in the sulcus are cut circularly with the aid of a perioto. Strong, long roots have to be loosened or luxated axially in both mesial and distal directions for a period of 30 seconds using a delicate lifting device / twister. In this process, care must be taken that no transversal movements are made. When extracting molars, the roots have to be separated before removal.

Indications for Use

The Bone Management® system Benex®-Control is indicated for the extraction of single-root tooth roots and tooth root fragments. In the case of multi-root tooth roots, these are to be separated before extraction.

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Insertion of the extraction screw

Depending on the case situation, the short or long extraction screw with the appropriate diameter is selected and can be carefully inserted into the tooth root using the jollying aid.

Note: Care must be taken that the jollying aid is pressed onto the extraction screw. This is the only way to ensure that there is a sufficient level of clamping force for safe transport into the oral cavity.

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Connection of the traction string with the protection cap

The images illustrate how the traction string is threaded into the protection cap. For ease of handling, it is an advantage to keep the traction string very short.

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Utilization

1. Fill the quadrant support with a fast-setting silicon or alginate and model a gap in the area of the buccal hole region.
2. Position the quadrant support with the gap over the root to be extracted, right-angle to the root axis onto the jaw.
3. Remove the quadrant support when the silicon or alginate starts setting. The entire setting process takes place extraorally.
4. After root luxation and the following insertion of the extraction screw, the repositioning of the quadrant support and the positioning of the extractor follow. The tooth root extraction is done according to the procedure in the point „Extraction of the tooth root“.

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Instructions for disassembly and assembly for the processing of the Benex extractor

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For optimal processing, the Benex Extractor must be completely disassembled prior to cleaning, disinfection and sterilization. To do this, first remove the knurled screw from the extractor by turning it counterclockwise off the thread.

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After that, the rotary segment can be removed from the main body by turning it clockwise or counterclockwise. A certain direction of rotation is not specified.

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By pulling or pushing the two cross-pins are to be detached from the body of the Benex extractor. It should be noted that the guide cover of the threaded pin has to be secured during this procedure to prevent the unwanted expansion of the compression spring.

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Now the guide cover can be removed carefully. The tensioning element should be fixed. After removing the compression spring, the tensioning element and the two O-rings can be removed.

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Subsequently, the Teflon segment plate can be released from the rotary segment by gently pulling.

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Now all individual parts of the extractor can be cleaned and disinfected according to valid processing instructions. After successful cleaning and disinfection and prior to packaging and sterilization of the Benex Extractor, the individual parts can be assembled. During assembly, it must be ensured that the O-rings have to be inserted prior to the cross pins in order to guarantee a perfect fit.

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When preparing the Benex-Control Extractor for sterilization, care must be taken to ensure that the ball does not engage in a depression of the segment plate but is in a semi-oblique position.

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Contraindications
Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

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Absolute contraindications
• Acute vascular diseases e.g. Myocardial infarction, severe heart failure, stroke
• Hematopoietic diseases such as acute leukemia and agranulocytosis
• Radiation and chemotherapy in the head and neck area
• Use in the case of bone fractures
• Dental ankylosis

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Relative contraindications
• Poor general condition of the patient (e.g. anemia)
• Acute infectious diseases as well as acute oral cavity infections
• Congenital tendency to infections (e.g. hypoimmunity, immunosuppressive therapy)
• Coagulopathy, anticoagulant therapy
• Poor condition and degree of damage of the adjacent teeth (e.g. loosening, unstable restorations)

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Special instructions for the extraction of curved tooth roots
• Unrooted teeth: The root must be well luxated intraalveolar. Very often, the curved root breaks during Benex extraction in the area of the most powerful curvature or at the level of the drill tip. In the second Benex drilling procedure or if the root curvature is very strong, the apical fragment can be extracted with the Benex extractor after the osteotomy according to the radius of the root curvature. The osteotomy in the apical root part is not necessary for the re-osteosynthesis of the alveolus or for the preservation of the alveolar ridge.
• Multi-rooted teeth: The root must be well luxated intraalveolar. It is recommended to extract the straight roots from the extractor first. This makes room for the conventional luxation of the curved root according to the radius of the curvature. The palatal roots of the maxillary molars are usually straight and can be removed with the extractor, while the buccal roots can be turned / pulled out with the lever. For removing the mandibular molars, the distal root is more suitable. The broad, often small mesial root can then be conventionally luxated / extracted in the distal direction. In most cases, this root is not suitable for Benex screw fixation.

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General instructions
Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

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