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Art.-No.: BS001 , BS002

DE

Application and safety instructions for the MEISINGER

Bone Scraper

Autologous Bone Collector

Art.-No.: BS001 , BS002

The use of autologous bone is the gold standard for bone augmentation. With the MEISINGER Bone Scrapers, autologous bone chips can be easily and safely collected and integrated into the present bone defect. Due to the immediate mixing with blood, the sterile gained bone material is ready for immediate use and can be placed on the defect. The curved shape of the bone chips creates a particularly large volume of augmentable material from bone and blood. Thanks to the completely detachable blade, the extracted bone can be easily removed and the individual parts can be safely processed after use. Another advantage: The exchangeable blades guarantee an optimal and constant cutting performance.

Indications

The Bone Management® Bone Scrapers are indicated for the manual harvest, collection and immediate application of autologous bone chips at the bone defect.

Utilization

Preparation

The two available versions of the Bone Scraper and the appropriate blades are made of surgical stainless steel. Before use, all individual parts have to be processed separately according to manufacturer information. Please note our instructions for processing. After processing, the appropriate blade is inserted into the selected Bone Scraper. Depending on the patient situation, the curved or straight version has to be selected.

Harvest and collection of the autologous bone

The Bone Scraper should be placed on the bone at a flat angle. Short, pulling movements with simultaneous medium pressure on the bone resulting in curved bone chips, which are collected in the sterile collection tray between the blade and the housing of the Bone Scraper. The bone chips obtained result in a high volume due to their curved shape, which significantly reduces the amount of augmentation material required and reduces the invasiveness of the intervention. The Bone Scrapers offer the possibility of obtaining both cortical and cancellous bone. The manual removal mixes the bone chips directly with patient blood and retains cellular components for optimal bone regeneration, such as osteocytes. The Bone Scrapers offer an alternative to conventional instruments such as trephines and bone mill. However, the autologous augmentation material obtained can also be used in combination with bone replacement materials.

Post-Processing

The rounded geometries and the complete dismantling allow for an uncomplicated cleaning of the Bone Scrapers, especially the collection tray. After use, all individual parts, in particular the blade used, are to be examined for signs of wear. As long as the items or the blades show no signs of wear, they can be reused. In case of doubt, however, the blades should be sorted out early to avoid damaging the bone. For the preparation of the individual parts including control, please observe our instructions for processing.

Note

Excessive pressure should be avoided when using the Bone Scraper. If sufficient cutting performance is no longer available, the used blade should be checked for wear and replaced if necessary.

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Contraindications

The use of the Bone Scrapers must comply with general medical and local, absolute and relative contraindications for dental surgery, which are relevant to the chosen augmentation procedure. In general: There must be no active local infections or pathological processes.


The absolute contraindications include severe internal diseases, diseases affecting bone metabolism, bleeding disorders, metabolic diseases, bone metabolism disorders (e.g. poorly controlled diabetes mellitus, osteoporosis), immunodeficiency and blood coagulation disorders.

The relative contraindications include psychiatric disorders, drug, alcohol and tobacco consumption, poor bone quality, poor oral hygiene of the patient, as well as head and neck radiotherapy, bisphosphonate therapy and other antiresorptive therapies.

NOTE: Contact of the collected bone with the saliva should be avoided. It is important to ensure that other structures such as soft tissue and nerves are not damaged.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

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