Additional Information on **Processing and Material Resistance of Autoclavable Anodized Aluminum Bur Blocks**

This anodized aluminum bur block requires careful treatment for processing.

- The Bur Blocks are suitable for manual cleaning and disinfecting as well as sterilization. Automated cleaning/disinfection with a disinfectant/WD (Washer-Disinfectant) should not be used due to material resistance.
- Remove instruments prior to cleaning and disinfecting and clean and disinfect in disassembled state.
- Cold-Sterile solution is not recommended by MEISINGER USA.

Note: If warm-air circulation is not evenly distributed, it may be helpful to place the block on its side. Discoloration and/or fading of the colored bur blocks occurring during sterilization, will not affect the quality of the block.

Warnings

The instructions for use described above, in particular for disinfection, cleaning and sterilization as well as speeds (if we add speed statement) must be observed. The products should only be used as intended. Failure to comply with safety instructions may result in damage and/or injury, such as heat necrosis, undesirable tissue dissection, nerve or tissue damage or infection

Information on material resistance:

Please consider the manufacturer's label for the cleaning and disinfecting agents: some are not suitable for aluminum.

When choosing the cleaning and disinfecting agents ensure that they do not contain the following ingredients:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 11, neutral/enzymatic, weak alkaline or alkaline cleaner recommended)
- organic solvents
- oxidizing agents
- halogens
- aromatic, halogenated hydrocarbons

In addition, please consider during selection of the detergents, that corrosion inhibitors, neutralizing agents, and/or rinse aids may cause potential critical remnants on the instruments. Please do not clean by use

of metal brushes or steel wool. Please do not expose to temperatures higher than 138 °C (280 °F)!

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing.

General Application and Safety Instructions for the Medical Device Produced by MEISINGER

MEISINGER medical products (dental, implantology, Bone Management[®], maxillary surgery, general surgery, chiropody devices) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.

It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the MEISINGER products for each individual case.

All Meisinger products have been developed for specific applications. Therefore, inappropriate use can result in

the premature wear and tear of the products and put patients and users at risk.

Disinfection, Cleaning and Sterilization

Unless there is explicit information to the contrary, all MEISINGER products are supplied in non-sterile packaging and, depending on the application, they need to be sterilized prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilized. Inappropriate cleaning and sterilizing of the instruments can result in the patient being infected with harmful bacteria.

You will find detailed instructions for the disinfecting, cleaning and sterilization of products in the Instructions for processing of medical devices produced by Hager & Meisinger GmbH. We would be happy to provide you with these instructions at your request. They are also available on the internet at www.meisinger.de.

Storage

During storage the products must be protected against dust, moisture and contamination. The maximum storage period must not be exceeded.

Products delivered non-sterile are marked with a standardized symbol "non sterile". If stored in accordance with regulations, the products can be stored for an unlimited period of time. For reasons of possible changes in material properties with regard to corrosion

resistance or the like, as well as the topicality of product design and suitability for archiving, the storage stability is set at a maximum of 10 years.

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