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Application and safety instructions for the MEISINGER

Crestal-Lift-Control Crestal-Lift-Control Basic

Crestal Approach Sinus Lift Systems

Art.-No.: CCLBA, BCL00

Crestal-Lift-Control is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially traumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Indications for Use

The Bone Management® systems Crestal-Lift-Control and Crestal-Lift-Control Basic are indicated for the performance of an internal sinus lift in case of a minimum vertical bone available about 5 mm. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Determination of the drilling depth / Selection of the stop sleeve

To determine the required drilling depth and selection of the appropriate stop sleeve, the analysis of the residual bone height is performed before the beginning of the surgery using an X-ray. The stop sleeve should be chosen so that about 2 mm residual bone remains to the maxillary sinus. At a residual bone height of e.g. 6 mm, the stop sleeve CL016 is selected to achieve a drilling depth of max. 4 mm.

Pilot drilling

The stop sleeve is placed onto the pilot bur and the initial drilling is performed. As cooling liquid, physiological, sterile saline solution should preferably be selected. We recommend to use the pilot bur intermittently, i.e. with up and down movements.


Extension drilling and access to the maxillary sinus

Afterwards the hole is extended with the aid of the crestal drills. For this purpose, the drills are used in ascending diameter until the required diameter for the implant insertion is reached. Until the penultimate extension drilling, the previously used stop sleeve (in the example CL016) is used. For the final extension drilling a 2 mm shorter stop sleeve (in the example CL014) is used for opening the access to the maxillary sinus. Due to the atraumatic construction of the crestal drills with the 4 cutting edges and the concave head, a conical bone flap can be shaped. This automatically causes the membrane to be raised by approx. 2 mm. Appropriate cooling has to be ensured during the entire drilling procedure. The crestal drills have to be used intermittently. In order to collect bone graft material, a speed of max. 50 rpm should be chosen.


Application of the depth gauge

The mobility of the bone flap can be checked with the depth gauge. Combined with the next smaller stop sleeve (in the example: CL013), the membrane can be lifted by another millimeter, if required.

Filling of the space with bone graft material

Afterwards, the new created space is filled up with bone graft material through the hole.

Spreading of the bone graft material

Using the bone graft condenser, the applied graft can be evenly spread within the sinus cavity. The application can be done by using the stop sleeve used for the final extension drilling.

Caution: The working part length of the bone graft condensers is 2.5 mm longer than the working part length of the crestal drills.

Implant insertion

If a one-stage procedure is indicated, the insertion of the selected implant can be carried out according to the manufacturer's instructions for use. The vertical minimum bone supply of approx. 5 mm allows for sufficient anchoring in the maxilla and thus also primary stability of the implant.

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Contraindications

The contraindications for the internal sinus lift are congruent with the contraindications for implantation. Therefore, relevant general conditions, diseases that may affect wound and bone healing, as well as a complete drug history should be performed. In addition, the following absolute contraindications should be noted:

- Residual bone height under 5 mm
- Not recommended in the case of complex morphology of the sinus floor (septum)
- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes

In case of smokers, care should be taken that there is no smoking one week before and at least 2 weeks after the procedure to support primary wound healing processes.

NOTE: A simple mucocoele in the preparation area (sinus floor) does not represent any contraindication and can usually be removed by means of aspiration.

CAUTION: Despite atraumatically designed instruments, care must be taken not to perforate the membrane.

ATTENTION: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

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