Absolute Kontraindikationen

Absolute contraindications

Contraindicaciones absolutas:

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- Insufficient supply of bone: Residual bone height < 12 mm / alveolar ridge width < 2 mm
- Implementation of the procedure directly after extraction or in the case of instable cortical bone
- Diseases that impair bone metabolism
- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes

NOTE:

Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® systems".

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² Largest working part diameter in 1/10 mm

Minimal diameter

External diameter

Pilotbohrer

Pilot burs

Fresas piloto

6 mm

8 mm

10 mm

12 mm

15 mm

Erweiterer

Expansion burs

Ensanchadores

L

Vorkörner

Initial bur

Fresa inicial

L

Diamantierte Säge

Diamond saw

Sierra diamantadas

L

Osteotomieinstrumente

Osteotomy instruments

Fresas para osteotomia

L

Handgriff

Handle

Mango

CBH0G

Fig.	A1001	D1001	A1003	D1003	186RF	231DC	859	859L	CBH0G
Shank ¹	206	206	206	206	204	204	204	314	
Size ²	010	018	023	030	018	070	018	012	
Length mm	15.0	15.0	15.0	15.0	12.0	0.3	10.0	12.0	
Minimal diameter			1.03	1.75					
External diameter			2.30	3.00					

¹ 204=RA, 206=RA XL, 314=FG ² Largest working part diameter in 1/10 mm Minimal diameter External diameter

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- Ridge widening of the mandible by selective spreading and angle modulation of the implant bed
- Controlled, minimally invasive procedures
- Maximal conservation of existing bone substance
- Suitable for preliminary implant cavity preparation for all current implant systems
- Optimal in combination with Split-Control

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Axial Angle Adjustment System

85FL005 - 0115

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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Horizontal-Control is an angle modulation system that enables easy and minimally invasive ridge widening, especially for the mandible with its approved special instruments. This is made possible by means of an angle modulation at the implantation axis determined by the supply of bone. Indications for the use of Horizontal-Control include "pear-shaped" atrophied or traumatized sections of the mandibular alveolar ridge that do not allow for a conventional implantation to be performed without appropriate Bone Management®. In this process, the cortical plate is straightened and shifted, condensing the bone tissue at the same time. For this purpose, the spreaders from the Bone Management® system, Split-Control, are required.

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First, an undersized pilot hole is made with the depth-marked pilot bur which determines the depth and direction of the implant bed (note diameter and depth setting).

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Using the conical expansion burs, the bone cavity is prepared for angle modulation.

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By carefully screwing in a spreader from one of the Bone Management® Split-Control systems, the implantation opening in the area of the cancellous bone is expanded and the adjacent bone tissue is condensed in a gentle manner. NOTE: Please also pay attention here to the application instructions for the spreaders.

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Depending on the accessibility within in the oral cavity, one of the cone cranks is screwed onto the handle and placed to form a perfect fit onto the spreader. By gradually changing the orientation of the instruments to the angle modulation of the implant axis in the direction of the bone deficit, the alveolar ridge can be widened by up to 8 mm.

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Depending on the indication, this process can be repeated with a spreader with a larger diameter so that the implant bed is properly prepared for inserting the implant.

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
The implant inserted immediately following this step holds the expanded gap open and stabilizes the cortical walls.

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First, an undersized pilot hole is made with the depth-marked pilot bur which determines the depth and direction of the implant bed (note

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If the existing conditions in the jaw do not

allow sufficient expansion, an internal vertical osteotomy can be performed in the area of the implant opening from distal and medial upwards (without making a flap). This creates a predetermined breaking point that breaks open the next time the opening is expanded. In this way, the cortical bone is shifted slightly outwards so the alveolar ridge can be spread to the required width.

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Using the wedge modulators, the moveable cortical plate is now bent carefully into the direction of the bone defect. The implant inserted immediately following this step holds the expanded gap open and stabilizes the cortical walls.

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