Application and safety instructions for the MEISINGER Lift-Control Systems

CONTENT

Language	Page
English (EN)	2-5
French (FR)	6-9
German (DE)	10-13
Italian (IT)	14-17
Spanish (ES)	18-21
Turkish (TR)	22-25

Application and safety instructions for MEISINGER Lift-Control Systems

Art.-No.: BTL00, CSL00, BLIPL

The Bone Management systems External-Lift-Control, Internal-Lift-Control and Lift-Control Plus offer instruments for the external sinus lift according to Boyne (direct sinus lift) as well as for the internal sinus lift according to Summers (indirect sinus lift).


INDICATIONS FOR USE

The Bone Management Systems are indicated for performing the external and internal sinus lift in cases of insufficient bone height in the maxillary posterior region. In the case of the internal sinus lift, augmentations of small vertical defects of 3 mm to 4 mm are possible with a residual bone height of ≥ 5 mm.

UTILIZATION


The specific system compositions as well as the recommended and maximum permitted speeds for the application are summarised in the article overview under „Content“. For proper, safe use, it is essential to observe these.

EXTERNAL SINUS LIFT (BTL00, BLIPL)


Lateral approach through the maxillary sinus wall

In the first step, the gingiva is folded down with a suitable incision in the area of the maxillary sinus so that the maxillary bone is exposed for further treatment. Using the spherical diamond and carbide instruments included in the set, a lateral window to the maxillary sinus is carefully prepared by thinning the bone in this area down to the membrane. The fine mucous membrane of the maxillary sinus (Schneiderian membrane) should not be perforated during the process.

Preparation of the Schneiderian membrane

The Schneiderian membrane is dissected off using the appropriate manual instruments. Depending on the anatomical conditions, the manual instruments are available in different sizes, shapes and angles to carefully loosen the membrane without injury and thus lift it cranially. The previously opened bone window of the maxillary sinus wall can now be carefully moved upwards into the maxillary sinus with the mucous membrane already detached.

Filling the cavity with bone substitute material


The cavity created in the maxillary sinus is now loosely filled with bone replacement material according to the manufacturer's instructions. The augmentation should hold the previously lifted bone window. The TL8 condenser can now be used to carefully compact the inserted augmentation. The still open lateral window is covered with the selected membrane according to the manufacturer's instructions. If necessary, the TL20 manual instrument can be used to position the membrane. The mucoperiosteal flap is then sutured in a saliva-proof and tension-free manner.

Implant insertion and further treatment


Depending on the patient's situation, implants can be placed in the same session (one-stage) or after the augmentation has healed (two-stage). The prerequisite for the single-stage procedure is the achievement of primary stability of the implant. This usually requires at least 3 mm of bone height. In the two-stage procedure, the implantation takes place only after complete ossification of the augmentation. After appropriate healing, the implant(s) is/are uncovered and the prosthetic treatment is placed. The corresponding manufacturer's instructions must be observed.

INTERNAL SINUS LIFT (CSL00, BLIPL)

Preparation of the implant site for indirect access to the maxillary sinus

The use of the pilot drilling can prevent the subsequently used twist drill from slipping, especially with a pointed alveolar ridge. A pilot hole is drilled using the depth-marked pilot drill B2001. The depth of this depends on the available vertical bone dimension, whereby the drilling should be carried out to a maximum of 1-2 mm below the sinus floor (X-ray monitoring). The hole is then extended according to the drilling protocol below.

Compaction of the jawbone and raising of the sinus floor

In the last step of the preparation, the bone is first compacted laterally and apically with the help of the corresponding elevator. The compacted bone is then driven cranially into the sinus cavity. At the same time, the Schneiderian membrane is raised, whereby care must be taken not to damage the membrane. Driving in the elevators should be done gently and by means of controlled hammer blows. This should be done by means of many short and gentle rather than less forceful hammer blows. If there is a strong cortical bone in the crestal area, it is recommended to widen it with the help of the pilot drill and thus facilitate the preparation with the elevators.

Depending on the individual patient situation, the choice of instruments used may vary.

The following applies:

Soft bone / condensation desired:

maximum drill diameter < maximum elevator diameter

Hard bone / condensation not desired:

maximum drill diameter > maximum elevator diameter

The diameter refers to the drill or elevator used in the previous preparation step.

In order to achieve an even higher compaction of the bone, if required, it is possible to prepare a smaller diameter hole and only achieve the desired diameter by using several elevators with increasing diameters. The preparation with the elevators is first carried out up to about 1 mm below the sinus floor. Only with the last elevator inserted is the sinus floor, including the bone compacted apically, driven cranially into the sinus cavity, with a simultaneous elevation of Schneiderian membrane.

Implant insertion

Possibility 1: Single-stage procedure

In the single-stage procedure, the newly created space is filled with the selected bone substitute material according to the manufacturer's instructions. Subsequently, the widening and fine preparation of the implant site can be carried out with the system's own instruments, as well as the insertion of the implant according to the manufacturer's instructions.

Possibility 2: Two-stage procedure

In the two-stage procedure, the first step is to fill the newly created space with the selected bone substitute material according to the manufacturer's instructions. The implantation is carried out according to the manufacturer's instructions in a second procedure as soon as the new bone formation is complete.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical measures must be observed. In addition, the contraindications for the implantation of an endosseous dental implant should be considered. The respective implant manufacturer provides information about these in the instructions for use of the selected implant system.

Absolute contraindications

- Complex morphology of the sinus floor (septum)
- Dentoalveolar growth not completed (exception: Cases where no dentoalveolar growth is expected, e.g. ectodermal dysplasia)
- Active infections and local pathological processes
- Taking medication that have an influence on bone healing (e.g. bisphosphonates)

COMMENT: A simple mucocele in the preparation area (sinus floor) is not a contraindication and can often be removed by suction.

Specific, absolute contraindications for the internal sinus lift

- Residual bone height < 5 mm

Relative contraindications

- Diseases that impair bone metabolism
- Drug and alcohol abuse
- Lack of patient cooperation
- Poor circulation
- Heavy work or active sports
- Mental state which may lead to non-compliance with the medical order
- Highly atrophied jaw

The removal of the Schneiderian membrane must be done very carefully, with low pressure. Perforation must be avoided. When using the elevators, take care not to change the axis of the periosteum, especially with very soft bone.

This product contains more than 0.1 % by mass of cobalt and therefore must be labelled as CMR substance class 1B (carcinogenic, mutagenic, reprotoxic) according to current regulations. It has been demonstrated that there is no increased risk of cancer or adverse reproductive or genetic effects when used as directed.

CAS: 7440-48-4
Contains hazardous substances:
Cobalt,
CAS: 7440-48-4

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disposal

When disposing of prosthetic components and instruments (at the end of their service life or after the specified shelf life has expired), care should be taken to ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national requirements (e.g. dual waste system).

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Manufacturer

Item number

Medical device

Batch number

Unique Device Identification

Do not re-use


Non-sterile

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Caution

Conformity mark

Do not use if the packaging is damaged

Manufacturing date

Packaging unit

CAS: 7440-48-4
Contains hazardous substances: Cobalt, CAS: 7440-48-4

Consult electronic instructions for use

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¹ 106=HR, 204=RA, 205=RA L

² Largest working part diameter in 1/10 mm

Minimal diameter

External diameter

CONTENT


Fig.	SL25	SL30	SL35	SL40	SL45	SL50	SL0	SL1
Shank ¹	-	-	-	-	-	-	-	-
Size ²	025	030	035	040	045	050	-	-
Length mm	16.0	16.0	16.0	16.0	16.0	16.0	100.0	133.0
	-	-	-	-	-	-	-	-
	2.5	3.0	3.5	4.0	4.5	5.0	-	-
Opt. speed rpm	-	-	-	-	-	-	-	-
Max. speed rpm	-	-	-	-	-	-	-	-
CE								
External-Life-Control	1x	1x	1x	1x	1x	1x	1x	1x
Internal-Life-Control								
Life-Control-Plus	1x	1x	1x	1x	1x	1x	1x	

¹ 106=HP, 204=RA, 205=RA L

² Largest working part diameter in 1/10 mm

Minimal diameter External diameter

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