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8FLB008 - 0421

Application and safety instructions for the MEISINGER
EN
Master-Pin-Control
Master-Pin-Control Basic

Pin System for Membrane Fixation

Art.-No.: BMP00, BMPBA


The Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5. Furthermore, they have a groove which increases the surface of the pins. In addition, due to the groove and a unique head-design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

Indications for Use

The Bone Management® systems Master-Pin-Control and Master-Pin-Control Basic are indicated for the fixation of apical membrane ends at the local bone in order to avoid micromobility of the membrane. They allow for the fixation of resorbable, non-resorbable and titanium-reinforced membranes.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Placing of bleeding points

For freshening the bone, bleeding points can be placed in the cortical bone in the region to be augmented using the included twist drills. If necessary, the bone can be pre-drilled with the aid of the initial drills. In case of hard bone, pre-drilling is recommended.

Application of the augmentation material

According to the respective manufacturer's instructions for use, the augmentation material is applied into the present defect and the selected membrane is adapted to the defect.

Picking up of the pins


The pins are located in a sterilizable box, which can be opened by a sliding mechanism on the side. With the aid of the pin holder, the pins are taken from the tray. For this, the pin holder is pressed vertically from above onto the pin heads. It is important to pay attention to a secure hold of the pins in the pin holder. The pins are not fixed in the box. Therefore, when opening the box as well as when picking up the pins with the pin holder, ensure vibration-free working. Note: Since the pins are implantable products, the pins should only be picked up using the pin holder to avoid touching them by hand.

Fixation of the membrane

To fix the membrane over the defect, the pins clamped in the pin holder are applied through the membrane into the cortical bone using a surgical hammer (e.g. HI070). By laterally tilting the pin holder, the holder can be released from the applied pin without negatively affecting the fixation of the pin in the bone. Vertical removal of the pin holder should be avoided as this may cause excessive pulling forces on the pin. In case of very soft bone, the use of the fixing holder is optionally indicated in order to fix the pins in the bone when releasing the pin holder. For this purpose, insert the uncotted working part of the fixing holder between the lamellae of the pin holder while disengaging the pin holder from the pin.

Important: Compression of the membrane must be avoided.

Removal of the pins

After a clinically confirmed, successful healing period, a small incision of the gingiva is made above the pin to be removed. Then, the pin can easily and very tissue-conservingly be loosened counterclockwise with a screwdriver; so that the loosened pin can be removed with the aid of tweezers or the fixing holder. When using the fixing holder care should be taken not to apply too much leverage, otherwise the flexible working part may be damaged.

NOTE: When removing the pin, note that there is no clamping effect provided between screwdriver and pin. The pins are a disposable product and are therefore not to be used again. The membrane should be handled according to the manufacturer's instructions for use.

Sterilization of the pins
The instruments and pins are supplied non-sterile and need to be cleaned, disinfected and sterilized prior to use.
For proper processing of the articles, please note the instructions for processing (cleaning, disinfection and sterilizing) of medical devices produced by Hager & Meisinger GmbH.

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Contraindications
Generally, medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Acute infection in or near the fixation site as well as pathological processes (e.g. symptoms like fever, local inflammation, abscess)
- Insufficient bone quality and/or quantity

Relative contraindications

- Known allergies and / or foreign body reaction to the implant material
- Systemic disorders
- Diseases that affect bone metabolism
- Medical treatments leading to progressive deterioration of bone (e.g. cortisone, immunosuppressive)
- Drug and alcohol abuse
- Lack of cooperation of the patient
- Poor circulation
- Hard physical labor or active sports
- Mental condition, which can lead to disregard of the medical order
- Highly atrophic jaw

IMPORTANT: Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection / dehiscence).

SEPARATE INSTRUCTIONS FOR PINS: Excessive physical activity and trauma affecting the membrane may result in premature failure of the pin due to loosening or breaking. The pins may only be used in accordance with the application and safety instructions. Any kind of modification in the procedure or on the product may impair the success of the treatment or damage the product. The pins are a disposable product and are therefore not be reprocessed or reused. Even with already used, but optically undamaged pins, the previous mechanical stress can lead to microscopic damage and thus to fatigue and produce failure. The pins have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the pins in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

NOTE ON MEMBRANE TECHNIQUE: Please note the information provided by the respective manufacturer of the membrane selected by you.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

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