Inhalt

Content Contenido

CE 0044

	Spiralbohrer für TCT09 Twist drills for TCT09 Taladros espirales para TCT09		Spiralbohrer für TCT10 Twist drills for TCT10 Taladros espirales para TCT10		Allport Allport Allport	Schraubendreher Screw driver Atornillador
Fig.	203RF	203RF	TC203	A2001	HM141F	31053
Shank ¹	205	205	204	204	205	205
Size	009	011	010	013	023	
Length mm	9.0	9.0	12.0	12.0		

¹204=RA, 205=RA L

²Largest working part diameter in 1/10 mm

Optional erhältlich: Schraubendreher
Optionally available: Screwdriver
Opción disponible: Atornillador

Art.-No. TL0T1

CE 0044

	Schrauben (je 7 St.) Screws (7 pcs. of each size) Tornillos (7 pzs.)		Schrauben (je 7 St.) Screws (7 pcs. of each size) Tornillos (7 pzs. c/u)	
Fig.	TCT09*	TCT09*	TCT09*	TCT09*
Length mm	7.0	10.0	13.0	16.0
	0.9	0.9	0.9	0.9
	1.1	1.1	1.1	1.1

ONLY in
Screw System TX Professional

CE 0044

	Schrauben (je 3 St.) Screws (3 pcs. of each size) Tornillos (3 pzs.)		Schrauben (je 3 St.) Screws (3 pcs. of each size) Tornillos (3 pzs.)	
Fig.	TCT10*	TCT10*	TCT10*	TCT10*
Length mm	7.0	10.0	13.0	16.0
	1.0	1.0	1.0	1.0
	1.3	1.3	1.3	1.3

 Minimaldurchmesser Minimal diameter Diámetro mínimo
 Außendurchmesser External diameter Diámetro externo

ONLY in
Screw System TX Professional

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Screw System TX Screw System TX Professional

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With the reuse of disposable products, the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

Hager & Meisinger GmbH | Hansemanstr. 10 | 41468 Neuss | Germany
Tel.: +49 (0) 21 31 20 12-0 | Fax: +49 (0) 21 31 20 12-222 | www.meisinger.de | info@meisinger.de

Meisinger USA, L.L.C. | 10200 E. Easter Avenue | Centennial, Colorado 80112 | USA
Tel.: +1 (303) 268-5400 | Toll free: +1 (866) 634-7464 | Fax: +1 (303) 268-5407
www.meisingerusa.com | info@meisingerusa.com

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
The osteosynthesis screws forming part of the Screw System TX sets are indicated for the universal fixation and mechanical stabilization of bone segments; for example, for the healing process in cases of vertical and horizontal augmentation grafts. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. In addition to the screws that are in the Screw System TX, the Screw System TX Professional contains screws in the lengths of 13.0 mm and 16.0 mm.

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
To optimize healing, the donor and recipient sites are prepared first with small dotted drillings (Bleeding Points). The bone transplant is adjusted to the bone site, which has been appropriately prepared, so that it will lie flat on the recipient bone.

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
The bone cylinder is drilled into by 0.9 mm using the twist drill (or 1.0 mm for rescue screws) and a drill hole is made up to the contralateral cortical bone of the recipient bone.

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
As a next step, only the bone cylinder is provided with a guiding hole, twist drill 1.1 mm (1.3 mm for rescue screw).

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Using the screw driver, a screw of the right length is used and the bone transplant is attached to the recipient bone. During this process, the screw should be gently tightened until the bone cylinder has been tightly fixed into place. Important: Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent tooth or tooth roots (risk of damage, infection / dehiscences). Depending on the case at hand, a decision should be made after the healing period as to whether the screws have to be removed or left in position. NOTE: The screws are a disposable product and are therefore not to be used again.

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NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

- Dentoalveolar growth which has not come to an end (exception: Cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes
- Cases in which the mechanical stability and the absence of mechanical movement of the transplant can not be achieved
- Insufficient bone supply (quality / quantity) and diseases that impair bone metabolism

NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

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