Application and safety instructions for the MEISINGER Split- and Crest-Control Systems

CONTENT

Language	Page
English (EN)	2-6
French (FR)	7-11
German (DE)	12-16
Italian (IT)	17-21
Spanish (ES)	22-26
Turkish (TR)	27-31

Application and safety instructions for the MEISINGER Split- and Crest-Control Systems

Art.-No.: CSPBA, CSP12, CSP15, BSPPR, CSPPL, CCR00

The present systems offer a comprehensive set of instruments for controlled ridge enlargement. They contain both special helical compaction and spreading instruments (spreaders) and / or flat horizontal spreaders that allow controlled and uniform spreading of horizontally resorbed bone. The helical spreading instruments (spreaders) also gently compact cancellous bone, while the horizontal spreaders allow the alveolar ridge to be spread by up to 5 mm. In this way, the jaw is optimally prepared for the insertion of an implant commonly available on the market.

INSTRUCTIONS FOR USE


The Bone Management Split-Control Systems are generally indicated for the expansion of the alveolar process with horizontal bone deficit (bone spreading) with simultaneous compaction of the bone if there is any low bone quality (bone condensing). This procedure is preferably used in the maxillary for density classes D3 and D4 (according to Misch).

The Crest-Control System is indicated for the expansion of the alveolar process by up to 5 mm with any horizontal bone deficit (bone splitting). This procedure is preferably used in the distal lower jaws.

UTILIZATION


The specific system compositions as well as the recommended and maximum permitted speeds for the application are summarised in the article overview under „Contents“. For proper, safe use, it is essential to observe these.

BONE SPREADING (CSPBA, CSP12, CSP15, CSPPL, BSPPR)


Setting the crestal relief incision

With the help of the diamond saw, a crestal relief incision is made, which supports the controlled expansion of the bone.

Setting the pilot bore

The pilot bore A1001 is placed along the planned implant axis. This is used to determine the depth and alignment of the implant bed. With a tapered alveolar ridge, the pilot drilling can be used in advance to prevent the subsequently used instruments from slipping off.

Use of the helical spreaders

A sequence of helical spreaders with increasing diameters is now screwed into the pilot bore so that the alveolar ridge is expanded along the previously created crestal relief incision and compacted at the same time (observe the sequence and depth marking of the instruments). The spreaders are gently screwed in by using the ratchet.

IMPORTANT: The helical spreaders are to be used with a torque of less than 35 Ncm. Furthermore, the minimum distance of 1.5 - 2 mm to neighbouring teeth and min. 3 mm to neighbouring implants must be observed.

For high bone density (D1-D2), the hole is prepared accordingly with the expanders before using the spreaders.

Insertion of the implants

The thus expanded implant bed now allows the insertion of a suitable implant according to the manufacturer's instructions. The significantly increased bone quality due to condensation also leads to optimised primary stability.

BONE CONDENSING (CSPBA, CSP12, CSP15, CSPPL, BSPPR)

In order to increase the density of the peri-implant bone before implant insertion, the bone condensing technique is used.

Setting the pilot bore

The pilot bore is placed along the planned implant axis. This is used to determine the depth and alignment of the implant bed. Depending on the planned implant diameter, use only the A1001 or the A1001 followed by the D1001.

With a tapered alveolar ridge, the pilot drilling can be used in advance to prevent the subsequently used instruments from slipping off.

Preparation for the use of the helical spreaders

The expanders now create an undersized, conical cavity that optimally prepares the bone for the use of the spreaders. Depending on the planned implant diameter, use only the A2003 or the A2003 followed by the D2003.

Use of the helical spreaders

The targeted use of the conical spreaders compacts the bone around the implant bed. For an implant bed to be compacted with a diameter of 2.7 mm or more, the A2003 expander is used together with the A2005 spreader. From a diameter of 3.3 mm, the D2003 expander is used together with the D2005 spreader. After the first spreader (A2005 or D2005) has been used, further spreaders of increasing diameter can be used without additional drilling. The spreaders are screwed in one after the other by using the ratchet.


Insertion of the implants

The thus compacted implant bed now allows the insertion of a suitable implant according to the manufacturer's instructions. The significantly increased bone quality due to condensation leads to optimised primary stability.


BONE SPLITTING (CCR00, CSPPL)

To achieve a widening of the alveolar ridge by up to 5 mm, the bone splitting technique is used.

Setting the crestal relief incision

With the help of the diamond saw, a crestal relief incision is made, which supports the controlled expansion of the bone.


Widening the crestal relief incision

The relief incision is widened with the help of the osteotomy instruments and thus prepared for the use of the horizontal spreaders. With a tapered alveolar ridge, the pilot drilling can be used in advance to prevent the subsequently used instruments from slipping off.


Inserting the horizontal spreaders

The closed horizontal spreaders are then positioned in the crestal relief incision.



Widening of the alveolar ridge

With the help of the SWOK I adapter, the horizontal spreaders are now spread slowly, alternately and equally, so that the alveolar ridge is widened evenly.

Insertion of the implants

The expanded alveolar ridge now allows the insertion of suitable implants according to the manufacturer's instructions. These can be inserted between the spreaders. Care must be taken not to damage the implant surfaces. Once the implant insertion is completed, the spreaders can be removed by closing them again with the help of the adapter. Cavities between the implants can, if necessary, be filled with bone substitute material according to the manufacturer's instructions.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered. In addition, the contraindications for the implantation of an endosseous dental implant should be considered. The respective implant manufacturer provides information about these in the user instructions of the selected implant system, among other things. In addition, the following contraindications must be observed:

Absolute contraindications

- Dentoalveolar growth not completed (exception: Cases where no dentoalveolar growth is expected, e.g. ectodermal dysplasia)
- Active infections and local pathological processes
- Application of the procedure directly after extraction or if the substantia compacta is not stable

Absolute contraindications for bone spreading / Bone Condensing

- Insufficient bone supply: Residual bone height < 12 mm / alveolar ridge width < 2 mm

Absolute contraindications for bone splitting

- Insufficient bone supply: Residual bone height < 10 mm / alveolar ridge width < 2 mm

Relative Contraindications

- Drug and alcohol abuse
- Lack of compliance by the patient
- Poor circulation
- Heavy labor or active sports
- Mental conditions that may cause patient's failure to adhere to physician's order
- Highly atrophic jaw

IMPORTANT: Care must be taken to protect the anatomical structures (safety distance min. 2 mm) and the course of the adjacent tooth/roots (risk of damage, infection, resorption).

The surgical and implantological procedure must always be carried out according to the current state of the art. The techniques described may only be used by trained professional users with sufficient experience with complex augmentation techniques. The following applies to all milling and drilling operations: To reduce the risk of bone overheating and thus necrosis, the instruments should be used with low contact pressure and constant cooling with sterile, physiological saline solution.

CAS: 7440-48-4

Contains hazardous substances:

Cobalt, CAS: 7440-48-4

This product contains more than 0.1 % by mass of cobalt and therefore must be labelled as CMR substance class 1B (carcinogenic, mutagenic, reprotoxic) according to current regulations. It has been demonstrated that there is no increased risk of cancer or adverse reproductive or genetic effects when used as directed.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national requirements (e.g. dual waste system).

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Manufacturer

Item number

Medical device

Batch number

Unique Device Identification

Do not re-use

Non-sterile

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Caution

Conformity mark

Do not use, if the packaging is damaged

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¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

Minimal diameter

External diameter

* **

CONTENT


Fig.	A2005	B2005	C2005	D2005	E2005	F2005	A I005	B I005	C I005	D I005	E I005	F I005
Shank ¹	-	-	-	-	-	-	-	-	-	-	-	-
Size ²	027	029	031	033	035	040	027	029	031	033	035	040
Lenght mm	12.0	12.0	12.0	12.0	12.0	12.0	15.0	15.0	15.0	15.0	15.0	15.0
	1.70	1.91	2.12	2.33	2.54	3.06	1.43	1.64	1.84	2.05	2.26	2.79
	2.7	2.9	3.1	3.3	3.5	4.0	2.7	2.9	3.1	3.3	3.5	4.0
Opt. speed rpm	-	-	-	-	-	-	-	-	-	-	-	-
Max. speed rpm	15	15	15	15	15	15	15	15	15	15	15	15
CE												
Crest Control												
Split Control Basic	1x			1x		1x						
Split Control 12 mm	1x	1x	1x	1x	1x	1x						
Split Control 15 mm							1x	1x	1x	1x	1x	1x
Split Control Plus	1x	1x	1x	1x	1x	1x						
Split Control Professional	1x	1x	1x	1x	1x	1x	1x	1x	1x	1x	1x	1x


Fig.	CA1RB	CA0RB	CA8RA	CA9RA	SW0L1	CARA4
Shank ¹	-	-	204	205	-	-
Size ²	-	-	-	-	-	-
Lenght mm	18.0	25.0	24.0	29.0	27.0	84.0
	-	-	-	-	-	-
	-	-	-	-	-	-
Opt. speed rpm	-	-	-	-	-	-
Max. speed rpm	-	-	15	15	-	-
CE						CE 0044
Crest Control					1x	1x
Split Control Basic		1x	1x			1x
Split Control 12 mm	1x	1x	1x			1x
Split Control 15 mm	1x	1x	1x			1x
Split Control Plus	1x		1x		1x	1x
Split Control Professional	1x	1x	1x	1x		1x

Fig.	VCD50
Width mm	6.0
Lenght mm	9.75
Height mm	10.6
Lifting Height mm	Max. 5.0
CE	
Crest Control	2x
Split Control Basic	
Split Control 12 mm	
Split Control 15 mm	
Split Control Plus	2x
Split Control Professional	

¹ 204=RA, 205=RA L, 206=RA XL ² Largest working part diameter in 1/10 mm Minimal diameter External diameter

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CE 0044

Crest Control	Ix		Ix	Ix	Ix	Ix							
Split Control Basic	Ix		Ix							Ix			
Split Control 12 mm	Ix	Ix		Ix	Ix	Ix					Ix		
Split Control 15 mm	Ix	Ix		Ix	Ix	Ix			Ix				Ix
Split Control Plus	Ix		Ix	Ix	Ix	Ix			Ix				Ix
Split Control Professional	Ix	Ix	Ix				Ix	Ix	Ix		Ix	Ix	Ix


Fig.	E 100I	F 100I	A 1003	D 1003	A2003	D2003	C3003	E3003	C3005	D3005	E3005	F3005	G2005	H2005
Gambo ¹	206	206	206	206	205	205	205	205	-	-	-	-	-	-
Taglia ²	020	025	023	030	023	030	028	032	031	033	035	040	045	050
Lungo mm	18.0	15.0	15.0	15.0	12.0	12.0	10.0	10.0	10.0	10.0	10.0	10.0	11.0	11.0
	-	-	1.03	1.75	1.30	2.03	2.00	2.42	2.30	2.51	2.72	3.24	3.65	4.17
	2.0	2.5	2.3	3.0	2.3	3.0	2.8	3.2	3.1	3.3	3.5	4.0	4.5	5.0
Velocidad óptima rpm	500-700	500-700	500-700	500-700	500-700	500-700	500-700	500-700	-	-	-	-	-	-
Velocidad máxima rpm	800	800	800	800	800	800	800	800	15	15	15	15	15	15

CE 0044

CE

Crest Control														
Split Control Basic														
Split Control 12 mm	Ix				Ix	Ix								
Split Control 15 mm			Ix	Ix										
Split Control Plus					Ix	Ix								
Split Control Professional	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix	Ix

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