Inhalt

Content Contenido

	Vorkörner Initial bur Fresa inicial	Diamantierte Sägen Diamond saws Sierras diamantadas		Osteotomieinstrumente Osteotomy instruments Fresas para osteotomía			Eindrehhilfen Jollying aids Soporte para enroscar				Ratsche Ratchet Carraca
CE 0044											
Fig.	186RF	231DC	231DC	HM331L	859	859L	CA1RB	CA0RB	CA8RA	CA9RA	
Shank¹	204	204	204	205	204	204	short	long	short	long	
Size²	018	070	100	010	018	010					
Length	12.0	0.3	0.3	5.5	10.0	12.0					

¹ 204=RA, 205=RA L ²Largest working part diameter in 1/10 mm

	Pilotbohrer Pilot burs Fresas piloto						Erweiterer Expansion burs Ensanchadores					
Fig.	A1001	B1001	C1001	D1001	E1001	F1001	A1003	D1003	A2003	D2003	C3003	E3003
Shank¹	206	206	206	206	206	206	206	206	205	205	206	206
Size²	010	013	015	018	020	025						
Length	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0	12.0	12.0	10.0	10.0
Minimal diameter							1.03	1.75	1.30	2.03	2.00	2.42
External diameter							2.30	3.00	2.30	3.00	2.80	3.20

¹ 206=RA XL ²Largest working part diameter in 1/10 mm Minimal diameter External diameter

	Spreader Spreaders Spreader						Erweiterer Expansion burs Ensanchadores						Spreader Spreaders Spreader					
Fig.	A1005	B1005	C1005	D1005	E1005	F1005	A2005	B2005	C2005	D2005	E2005	F2005	C3005	D3005	E3005	F3005	G2005	H2005
Length	15.0	15.0	15.0	15.0	15.0	15.0	12.0	12.0	12.0	12.0	12.0	12.0	10.0	10.0	10.0	10.0	11.0	11.0
Minimal diameter	1.43	1.64	1.84	2.05	2.26	2.79	1.70	1.91	2.12	2.33	2.54	3.06	2.30	2.51	2.72	3.24	3.65	4.17
External diameter	2.70	2.90	3.10	3.30	3.50	4.00	2.70	2.90	3.10	3.30	3.50	4.00	3.10	3.30	3.50	4.00	4.50	5.00

Minimal diameter External diameter

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Split-Control Professional

Bone Spreading and Condensing System

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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
The Split-Control Professional system forming part of the Bone Management® Series makes innovative instruments available specifically for bone spreading and bone condensing. The system is indicated for the widening of the alveolar process in the case of a horizontal bone deficit with simultaneous lateral bone condensing, which provides for better primary stability for the implantation. The procedure is preferred for use in the maxilla at density classes D3 and D4 (according to Misch). Bone spreading is a minimally invasive technique for the dilation of horizontally resorbed bone in a gentle and controlled manner. By means of targeted bone condensing, the bone is additionally condensed and the primary stability of the inserted implant is increased. By using the screw-shaped spreaders (spreading and condensing instruments), the bone is thereby optimally prepared for the insertion of an implant. Split-Control Professional contains an extensive set of instruments and a wide selection of spreaders in order to be able to selectively deal with a variety of specific indications.

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
Using a diamond saw, a crestal relief cut is made that promotes the controlled dilation of the bone.

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
The initial bur is used to prevent slippage of the pilot bur that is subsequently used.

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The successive deployment of a series of spreaders with increasing diameters leads to bone spreading and condensing (order of diameter and codification of the instruments must be respected). With the help of the ratchet the spreaders may be rotated in cautiously. Torque under 35 Ncm is recommended. Furthermore, a minimum distance of 1.5 - 2 mm away from adjacent teeth and 3 mm away from adjacent implants should be adhered to. In the case of high bone density (D1-D2), drilling takes place before using the spreaders with expansion burs that have been prepared accordingly.

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
By having widened the implant cavity in this manner it is now possible to place the suitable implant. Furthermore bone rigidity, achieved by bone condensation, may lead to greater primary stability.

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In the case of a narrow alveolar ridge, the initial bur can be used to prevent slippage of the pilot bur that is subsequently used. With the aid of the pilot bur, sub-dimensional pre-drilling can be performed afterwards. Please pay attention to the drilling depth and diameter when selecting the related instruments.

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A sub-dimensional bone cavity for bone densities D3/D4 is then created using the conical expansion burs.

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For bone condensing purposes of implant cavities of more than 2.7 mm diameter, use the expansion bur A1 and the spreader A1. For implant cavities exceeding 3.3 mm diameter, the expansion bur D1 and the spreader D1 must be used. After placement of the primary spreader (A1 or D1), the following spreaders in ascending order can be used without further drilling with pilot drills. In non-atrophied bone, the use of the spreaders leads to condensation of the bone and horizontally resorbed areas of the alveolar ridge are expanded in addition. With the help of the ratchet, the spreaders are rotated in.

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The condensed and spread bone may now provide greater primary stability once an implant has been placed.

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• Insufficient supply of bone: Residual bone height < 12 mm / alveolar ridge width < 2 mm
• Implementation of the procedure directly after extraction or in the case of instable cortical bone
• Diseases that impair bone metabolism
• Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
• Active infections as well as pathological processes
NOTE: Related contraindications and further instructions deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

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