## Inhalt Split-Control 12 mm / Split-Control 15 mm

Content Split Control 12 mm / Split Control 15 mm Contenido Split Control 12 mm / Split Control 15 mm

	Vorkörner Initial bur Fresa inicial	Diamantierte Sägen Diamond saws Sierras diamantadas		Kronentrenner Crown cutter Fresa para separar coronas	Osteotomieinstrumente Osteotomy instruments Fresas para osteotomia		Ratsche Ratchet Carraca	Eindrehhilfen Jollying aids Soporte para enroscar		
Fig.	186RF	231DC	231DC	HM31C	HM33IL	859	859L	CARA4	CA1RB	CA0RB
Shank <sup>1</sup>	204	204	204	204	205	204	204			
Size <sup>2</sup>	018	070	130	012	010	018	010			
Length mm	12.0	0.3	0.3	4.1	5.5	10.0	12.0	84.0		

<sup>1</sup> 204=RA, 205=RA L <sup>2</sup> Largest working part diameter in 1/10 mm

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<sup>1</sup> 205=RA L, 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

	Spreader Spreaders Spreader					
Fig.	A2005	B2005	C2005	D2005	E2005	F2005
Size	027	029	031	033	035	040
Length mm	12.0	12.0	12.0	12.0	12.0	12.0
▲	1.70	1.91	2.12	2.33	2.54	3.06
▲	2.70	2.90	3.10	3.30	3.50	4.00

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<sup>1</sup> 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

	Spreader Spreaders Spreader					
Fig.	A1005	B1005	C1005	D1005	E1005	F1005
Size	027	029	031	033	035	040
Length mm	15.0	15.0	15.0	15.0	15.0	15.0
▲	1.43	1.64	1.84	2.05	2.26	2.79
▲	2.70	2.90	3.10	3.30	3.50	4.00

▲ Minimal diameter      ▲ External diameter

## Inhalt Split-Control Basic

Content Split Control Basic Contenido Split Control Basic

	Vorkörner Initial bur Fresa inicial	Diamantierte Säge Diamond coated saw Sierras diamantadas	Ratsche Ratchet Carraca
Fig.	186RF	231DC	CARA4
Shank <sup>1</sup>	204	204	
Size <sup>2</sup>	018	100	
Length mm	12.0	0.3	84.0

<sup>1</sup> 204=RA <sup>2</sup> Largest working part diameter in 1/10 mm

	Pilotbohrer Pilot burs Fresas piloto	Spreader Spreaders Spreader			Eindrehhilfe Jollying aid Soporte para enroscar
Fig.	A2001	A2005	D2005	F2005	CA0RB
Shank <sup>1</sup>	205				
Size <sup>2</sup>	013	027	033	040	long
Length mm	12.0	12.0	12.0	12.0	
▲		1.70	2.33	3.06	
▲		2.70	3.30	4.00	

<sup>1</sup> 205=RA L <sup>2</sup> Largest working part diameter in 1/10 mm

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- Controlled spreading of horizontally resorbed bone
- Gentle bone condensing by means of non-cutting Archimedes screws
- High level of primary stability of the inserted implants
- Effective in increasing the width of narrow alveolar ridges
- Preparation of the bone for all standard forms of implant systems

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## Split-Control Basic 12 mm Split-Control 12 mm Split-Control 15 mm

Bone Spreading and Condensing System

85FLBM10 - 0315

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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Meisinger since 1888

# Split-Control Basic 12 mm Split-Control 12 mm Split-Control 15 mm

## Bone Spreading and Condensing System

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
The Split-Control systems forming part of the Bone Management® Series make innovative instruments available specifically for bone spreading and bone condensing. The systems are indicated for the widening of the alveolar process in the case of a horizontal bone deficit with simultaneous lateral bone condensing which provides for better primary stability for the implantation. The procedure presented is preferred for use in the maxilla at density class D3 and D4 (according to Misch). Bone spreading is a minimally invasive technique for the dilation of horizontally resorbed bone in a gentle and controlled manner. By means of targeted bone condensing, the bone is additionally condensed and the primary stability of the inserted implant is increased. By using the screw-shaped spreaders (spreading and condensing instruments), the bone is optimally prepared for the insertion of an implant.

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
Using a diamond saw, a crestal relief cut is made that promotes the controlled dilation of the bone.

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
The pilot bur is used at the areas where the implants should later be inserted. In order to avoid slippage, the initial bur can be primarily used.

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
The successive deployment of a series of spreaders with increasing diameters leads to bone spreading and condensing (order of diameter and codification of the instruments must be respected). With the help of the ratchet the spreaders may be rotated in cautiously. Torque under 35 Ncm is recommended. Furthermore, a minimum distance of 1,5 - 2 mm away from adjacent teeth and 3 mm away from adjacent implants should be adhered to. In the case of high bone density (D1-D2) drilling takes place before using the spreaders with expansion burs that have been prepared accordingly.

The successive deployment of a series of spreaders with increasing diameters leads to bone spreading and condensing (order of diameter and codification of the instruments must be respected). With the help of the ratchet the spreaders may be rotated in cautiously. Torque under 35 Ncm is recommended. Furthermore, a minimum distance of 1,5 - 2 mm away from adjacent teeth and 3 mm away from adjacent implants should be adhered to. In the case of high bone density (D1-D2) drilling takes place before using the spreaders with expansion burs that have been prepared accordingly.

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
By having widened the implant cavity in this manner it is now possible to place the suitable implant. Furthermore bone rigidity, achieved by bone condensation, may lead to greater primary stability.

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
The initial bur can be used to prevent slippage of the pilot bur that is subsequently used. With the aid of the pilot bur, sub-dimensional pre-drilling can be performed afterwards. Please pay attention to the drilling depth and diameter when selecting the related instruments.

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
A sub-dimensional bone cavity for bone densities D3/D4 is then created using the conical expansion burs.

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
For bone condensing purposes of implant cavities of more than 2.7 mm diameter, use the expansion bur A1 and the spreader A1. For implant cavities exceeding 3.3 mm diameter, the expansion bur D1 and the spreader D1 must be used. After placement of the primary spreader (A1 or D1), the following spreaders in ascending order can be used without further drilling with pilot drills. In non-atrophied bone, the use of the spreaders leads to condensation of the bone and horizontally resorbed areas of the alveolar ridge are expanded in addition. With the help of the ratchet, the spreaders are rotated in.

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The condensed and spread bone may now provide greater primary stability once an implant has been placed.

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- Insufficient supply of bone: Residual bone height < 12 mm / alveolar ridge width < 2 mm
- Implementation of the procedure directly after extraction or in the case of instable cortical bone
- Diseases that impair bone metabolism
- Dentoalveolar growth which has not come to an end (exception: Cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes

NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

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