Inhalt Transfer-Control

Content Contenido

	Osteotomieinstrumente Osteotomy instruments Fresas para osteotomía			Diamantierte Säge Diamond saw Sierra diamantada
Fig.	HM33IL	859L	859	231DC
Shank ¹	205	204	204	204
Size	010	010	018	070
Length mm	5.5	12.0	10.0	0.3

¹ 204=RA, 205=RA L ²Largest working part diameter in 1/10 mm

	Trepane Trephines Trépanos			Radfräser Wheel cutters Fresas rueda			Lagerfräser Ablative burs Fresas ablativas		
Fig.	T229L	T229L	T229L	TC084	TC084	TC084	TC21X	TC21X	TC21X
Shank ¹	104 / 205	104 / 205	104 / 205	204	204	204	104 / 205	104 / 205	104/205
Size ²	040	050	060	004	005	006	040	050	060
Length mm	14.0	14.0	14.0	4.0	5.0	6.0	14.0	14.0	14.0
	4.0	5.0	6.0						
	5.0	6.0	7.0	10.0	10.0	10.0			

¹ 104=HP, 204=RA, 205=RA L ²Largest working part diameter in 1/10 mm External diameter Internal diameter

Inhalt Transfer-Control Basic

Content Contenido

	Diamantierte Kugel Diamond coated ball Bola diamantada	Diamantierte Säge Diamond saw Sierra diamantada
Fig.	801	231DC
Shank ¹	205	204
Size	033	070
Length mm		0.3

¹ 204=RA ²Largest working part diameter in 1/10 mm

	Trepane Trephines Trépanos		Lagerfräser Ablative burs Fresas ablativas	
Fig.	T229L	T229L	TC21X	TC21X
Shank ¹	205	205	205	205
Size ²	050	070	050	070
Length mm	14.0	14.0	14.0	14.0
	5.0	7.0		
	6.0	8.0		

¹ 205=RA L ²Largest working part diameter in 1/10 mm External diameter Internal diameter

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- Easy and safe transplantation of bone cylinders
- Instruments that are intelligently matched to one another
- Perfectly fitting bone cylinders
- Precise and predictable application due to depth marks
- Optimally suitable for vertical and horizontal augmentative grafts
- Rapid vitalization and tissue repair

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Transfer-Control Transfer-Control Basic

Bone Replacing System

85FL002 - 0115

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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The MEISINGER Transfer-Control systems are indicated for the intraoral extraction of bone cylinders as well as for horizontal and vertical bone deficits in the maxilla and the mandible. With these systems, the most important conditions for the healing of horizontal and vertical augmentative grafts, and a congruent and fresh bone recipient site can be created in an easy and controlled manner. The tools are matched to each other in such a way that the outer diameters of the ablative burs and the diameters of the wheel cutters correspond to the inner diameter of the trephines. In this way, for the most part, a tug back is aimed at when inserting the bone cylinder that can be reinforced with osteosynthesis screws, which are separately available. Perfectly transplanted bone cylinders of this type result in a bone site that is capable of implantation after only approx. 3–4 months due to quicker vitalization and healing.

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
In cases of severe horizontal alveolar ridge resorption (where even bone-spreading techniques cannot be deployed) **horizontal augmentative grafts** are required for compensation purposes.

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Depending on the implant diameter required, initially a bone cylinder of the correct diameter is extracted using an appropriate trephine, e.g. from the retromolar region of the mandible. If necessary, the osteotomy disk or diamond can be used to create a horizontal cut allowing the bone cylinder to be taken from the donor site.

NOTE: Should no initial bur be used, the trephine should be initially used rotating toward the left (fewer rotations).

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Next, using the ablative bur corresponding to the trephine used previously, a bone site is ground on the horizontally resorbed alveolar process at the tooth axis. The surgical technique and the properties of the ablative bone bur lead to a hyperaemic environment, which accelerates the wound healing in post transplantation. The donor site may be filled with collagen or other grafting material. The intervention is terminated by a current gingival suture.

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The extracted bone cylinder is placed into the recipient bone site to compensate the horizontal bone deficit.

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
Depending on the case situation at hand, the bone cylinder can be attached with a fixation screw from the Screw System TX (Art.-No. BTX00, BTXPR),

which optimally compliments Transfer-Control.

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
In cases of severe vertical alveolar ridge resorption, **vertical augmentative grafts** are required for compensation. In this case the appropriate wheel-cutter creates a congruent cavity in the crestal alveolar ridge. One should take note of the correlating trephines to be used. Prepare an autogenous bone graft in the same manner as described for horizontal defects.

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