**General application and safety instructions for the
medical device produced by MEISINGER**

CONTENT

Language	Page
English (EN)	2
French (FR)	3
German (DE)	4
Italian (IT)	5
Spanish (ES)	6

General application and safety instructions for the medical device produced by MEISINGER

- MEISINGER medical products (dental, implantology, Bone Management®, maxillary surgery, general surgery, chiropody devices) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.
- It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the MEISINGER products for each individual case.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- For information on possible applications of MEISINGER products, please refer to the manufacturer's specifications (see catalogue)
- Courses providing further training on the possible uses of MEISINGER products are organised on a regular basis. You will find relevant information about these events as well as further information on Meisinger products on the internet at www.meisinger.de.
- All Meisinger products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

Warnings

The described instruction for use, in particular for pressure force, cooling, disinfection, cleaning and sterilisation as well as speeds, must be observed. The instruments should only be used as intended. Failure to comply with the safety instructions may result in damage to the drive and / or injury, such as heat necrosis, undesirable tissue dissection, nerve or tissue damage or infection.

Application

- Care must be taken to use only functional, technically and hygienically flawless, maintained and cleaned turbines, straight and contra-angle handpieces.
- Rotating instruments must be clamped as low as possible and must be checked for secure seating of the product in the turbine, handpiece or angled piece before use.
- The rotating instruments may only be used with transmission instruments (e.g. handpieces, contra-angles) that are approved as medical devices and suitable for the intended purpose. The transfer instruments must have a standardised interface for rotating instruments. The manufacturer's application and safety instructions must be observed.
- The instruments must be rotating when applied on material. They should not be placed on material and then brought to rotation.
- Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- Inappropriate use of the products leads to increased risk and bad work results.
- When working with dry materials, it is recommended to use a suction cleaning device.
- In particular, users of hand tools should take care to use them gently and with consideration.
- The instruments shall only be touched with protective gloves.
- Thermal bone damage caused by rotating tools (e.g. pilot burr, conical burr, expansion burr) must at all times be avoided (user training, working at low speed and with sufficient cooling).
- During intraoral application attention has to be made to the fact that the products are protected against aspiration or falling on the floor.

Contraindications

The products must not be used outside their intended purpose and range of indications. The specified maximum speeds must not be exceeded.

Use of pressure

- Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- During polishing, excess pressure can lead to heat generation.
- Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

Cooling

- In order to avoid excessive heat generation during preparation, a sterile water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

Disinfection, cleaning and sterilization

- Unless there is explicit information to the contrary, all MEISINGER products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Improper cleaning and sterilization of instruments exposes the patient to an increased risk of infection.
- You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the processing of medical devices produced by Hager & Meisinger GmbH. We would be happy to provide you with these instructions at your request. They are also available on the internet at www.meisinger.de.

Speed recommendations for rotary instruments

- Following the instrument-specific speed recommendations produces the best results.
- Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.
- Please consult the manufacturer's information (see catalogue, speed overview or www.meisinger.de) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.

• Generally, the following rules apply:

- The larger the working part of an instrument the lower the speed
- The larger the working part of an instrument the greater the force of pressure
- Instruments with a maximum speed of 300,000-450,000 rpm are suitable for micro-motor hand pieces and turbines with stable ball bearings. Not to be recommended for turbines with air bearings.
- Instruments with a maximum speed of 30,000 to 160,000 rpm are suitable for micro-motor hand pieces or technical hand pieces up to the indicated rpm. Not to be recommended for turbines.
- Surgical instruments: suitable for geared down micro-motor hand and angle pieces 10:1 with stable ball bearings. The instruments are used with physiological, possibly with sterile external cooling or internal cooling if an appropriate handpiece is used.

Discarding worn instruments and parts

- MEISINGER products can principally be reused several times – unless specifically indicated and labelled otherwise. Rotating instruments are subject to wear. The option of and accountability for multiple use of a product and the frequency of application is solely the decision and own responsibility of the treating clinician based on the application in each case and the possible wear of the products. If in doubt, the products should always be sorted out early and replaced.
- Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- Instruments that are bent and/or do not run true should be discarded forthwith.
- With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

Storage

During storage the products must be protected against dust, moisture and contamination. The maximum storage period must not be exceeded.

- Products delivered in sterile condition (e.g. implants, singles) have a limited shelf life. The shelf life of the products is indicated by the standardised symbol „hourglass“ and must be observed. If the products are used after the shelf life has expired or if the packaging is damaged, sterility cannot be guaranteed.
- Products delivered non-sterile are marked with a standardised symbol „non sterile“. If stored in accordance with regulations, the products can be stored for an unlimited period of time. For reasons of possible changes in material properties with regard to corrosion resistance or the like, as well as the topicality of product design and suitability for archiving, the storage stability is set at a maximum of 10 years.
- Non-sterile products with a minimum shelf life (e.g. polishers) have a limited shelf life. The shelf life is indicated by the standardized "hourglass" symbol and must be observed. If the products are used after the shelf life has expired, the material durability cannot be guaranteed due to possible changes in the material properties.

Additional instructions regarding the use of trepans

- When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

Additional instructions for the use of diamond articles

Diamond articles can be used several times in the dental field with appropriate processing. Used in surgery, i.e. surgically invasive application to the bone, the diamond instruments are to be used as disposable products.

Additional information for the implementation of the SINGLES (supplied sterile single use product)

- SINGLES are delivered in sterile form (sterilization with gamma rays, radiation intensity at least 25kGy). The inner packaging consists of a special tubular bag and creates a sterile barrier of the products. The tubular bag is in a folding box, which is provided with a label and has the marking of the product. The products are intended for single use.
- If the inner packaging is damaged and/or unsealed, the sterility of the product cannot be guaranteed, resulting in an increased risk of contamination. In such a case, the product is not allowed to be sterilized and/or used and must be disposed. After the expiry of the indicated shelf life of the sterilization (hourglass symbol), the products may also not be used.
- To remove the instruments, tear open the pouch at the side perforation. The instruments should not be pushed out through the sterile packaging.

Additional instructions for the use of rotating root canal instruments

Root canal instruments are only to be used for preparing and shaping the coronal part of the root canals. The mechanical extension takes place in the coronal third of the root canal and the root canal entrance. The instruments are used for the mechanical preparation of the access to the root canal entrances by completely removing the pulp wax. The use of the instruments is only permitted in the straight portion of the root canals.

Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. Dual Waste System).

Further comments:

- Due to legal requirements, returned goods can generally only be accepted in a sealed original packaging with the indication of the complete lot number. This can be found on the product packaging.

Hager & Meisinger GmbH
Hansemannstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

(*)
This product contains more than 0.1 % by mass of cobalt and therefore must be labelled as CMR substance class 1B (carcinogenic, mutagenic, reprotoxic) according to current regulations.
It has been demonstrated that there is no increased risk of cancer or adverse reproductive or genetic effects when used as directed.

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