



Universal Extraoral Polishing & Shaping Kit

according to Dentsply Sirona®

NEW!



Item No. #DSU21
Universal Extraoral Polishing & Shaping Kit
according to Dentsply Sirona®

At a glance

- Comprehensive set of burs and polishers for all your shaping and polishing needs
- Contour and polish in the pre-sintered and post sintered state
- Create stunning restorations
- Designed for all zirconia and glass ceramic materials

More Information

In conjunction with Dentsply Sirona®, MEISINGER has created this Universal Extraoral Polishing and Shaping Kit to provide dental professionals with an excellent, simple, and versatile solution for all of their modern CAD/CAM contouring and polishing needs.

Included are MEISINGER's best-selling twist polishers, diamonds, and carbides to contour and polish restorations with a unique, personal flare. This kit is ideal for all CAD/CAM restorative materials.

Materials

Intended for use with all CAD/CAM materials including:

- CEREC Tessera™ (registered trademark of DENTSPLY SIRONA)
- CEREC MTL™ (registered trademark of DENTSPLY SIRONA)


MEISINGER USA

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
customerservice@meisingerusa.com • www.meisingerusa.com

@meisingerusa

Meisinger
since
1888

GLASS CERAMIC



CONVENIENT ENGRAVINGS

The Universal Extraoral Polishing and Shaping Kit provides clearly laser-engraved sections for instruments to be used on Glass Ceramics, Zirconia, and items that can be used on both materials.

ZIRCONIA

Items that overlap the Glass Ceramic and Zirconia sections can be used on both materials.

Material(s) Legend:

■ Glass Ceramic
 ■ Glass Ceramic & Zirconia
 ■ Zirconia

#DSU21 - KIT CONTENTS						
	Reorder #	Usage	Diameter (1/10mm)	Length (mm)	Shank	Pcs/Kit
	936SC-200-HP	Separation	200	0.2	HP	1
	9736G-150-HP	Adjusting	150	3.0	HP	1
	850F-025-HP	Adjusting	025	8.0	HP	1
	9752G-170-HP	Pre-Polishing - Axial Walls	170	2.5	HP	1
	9681P-220-HP	High Shine - Axial Walls	220	3.0	HP	1
	805S-018-HP	Characterization	018	2.5	HP	1
	9771G-170-HP	Pre-Polishing	170	1.6	HP	1
	9771M-170-HP	Polishing	170	1.6	HP	1
	9771F-170-HP	High Shine Polishing	170	1.6	HP	1
	9769F-170-HP	Final Polish	170	1.6	HP	1

#DSU21 - KIT CONTENTS (continued)						
	Reorder #	Usage	Diameter (1/10mm)	Length (mm)	Shank	Pcs/Kit
	HM488FX-016-HP	Separation / Characterization	016	–	HP	1
	HM79EX-040-HP	Smooth Sprue	040	–	HP	1
	9617V-220-HP	Adjusting	220	–	HP	1
	DCA06-170-HP	Pre-Polish - Axial Walls	170	2.5	HP	1
	DCA12-170-HP	High Shine - Axial Walls	170	2.5	HP	1
	DCA04-040-HP	Pre-Polish - Occlusal Surface	040	–	HP	1
	DCA10-040-HP	High Shine - Occlusal Surface	040	–	HP	1

COMPLETE SYSTEMS - see page 08-09 for more information

	Complete System (Kit, Organizer, & Refills)	Complete System Refill (Organizer & Refills)
	#DSU21A	#DSU21AR

MEISINGER USA

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 customerservice@meisingerusa.com • www.meisingerusa.com

@meisingerusa



Visit www.meisingerusa.com to browse and shop our entire catalog of products.