

New
Luster[®]

Lithium Disilicate Polishers by MEISINGER

Ideal for Today's New Dental Restorations

Luster[®] for Lithium Disilicate

Extraoral Adjusting and Polishing Kit

This kit features the MEISINGER specially designed polishers for Lithium Disilicate restorations. The 9736G wheel allows for adjustments, and the two-step polishing system will re-glaze the restoration.

Adjusting and Polishing Kit # LUS80



ALL FOUR KITS

rated "CR Choice"
by Clinicians Report[®] for
best polishing on
Lithium Disilicate
and Zirconia.

Luster[®] for Lithium Disilicate

Intraoral Adjusting and Polishing Kit

This kit contains specially formulated diamonds to make adjustments to Lithium Disilicate as well as the two-step polishing system to re-glaze Lithium Disilicate intraorally after the restoration has been placed.

Intraoral Adjusting Kit # LUS100



* "An independent, non-profit, dental education and product testing foundation, Clinicians Report[®], February 2013. For the full report to be sent directly to you, please contact MEISINGER USA customer service."

MEISINGER USA, L.L.C.

10200 E. Easter Ave. • Centennial • Colorado • 80112 • USA

Tel.: +1 (303) 268-5400 • Fax: +1 (303) 268-5407

E-Mail: info@meisingerusa.com • <http://meisingerusa.com> • www.bone-management.com • www.occlusalrouter.com

Meisinger
125
YEARS

Luster®

Zirconia Polishers by MEISINGER

Ideal for Today's New Dental Restorations

Luster® for Zirconia

Extraoral Adjusting and Polishing Kit

This kit was designed to adjust and re-glaze monolithic zirconia restorations such as BruxZir® crowns.

Adjusting and Polishing Kit # LUS85



ALL FOUR KITS

rated "CR Choice"
by Clinicians Report® for
best polishing on
Lithium Disilicate
and Zirconia.

Luster® for Zirconia

Intraoral Adjusting and Polishing Kit

This kit, with specially formulated polishers, was designed to adjust and re-glaze monolithic zirconia restorations such as BruxZir® crowns intraorally.

Intraoral Adjusting Kit # LUS91



* "An independent, non-profit, dental education and product testing foundation, Clinicians Report®, February 2013. For the full report to be sent directly to you, please contact MEISINGER USA customer service."

MEISINGER USA, L.L.C.

10200 E. Easter Ave. • Centennial • Colorado • 80112 • USA

Tel.: +1 (303) 268-5400 • Fax: +1 (303) 268-5407

E-Mail: info@meisingerusa.com • <http://meisingerusa.com> • www.bone-management.com • www.occlusalrouter.com

Meisinger
125 YEARS

BruxZir® is a registered trademark of Glidewell Laboratories.