

CF-Diamonds

For easy and precise preparation

The MEISINGER CF-Diamonds have a unique diamond finish with two different grain sizes. The tip of the instrument has a fine grain, the tapered area of the working part is covered with a coarse diamond grain. The instrument is particularly suitable for the preparation of anterior teeth. With the finely coated tip, a perfect preparation margin can be defined without any change of the instrument.



At a glance

- Suitable for preparation in the anterior region
- Perfect definition of the preparation margin
- Saved material removal to the incisal edge of the tooth
- Gentle preparation of the tooth
- Long service life

Our sales force would be pleased to
answer your questions: +49 2131 2012-303

Meisinger
since
1888

