

Composite Singles

Single-use polishers to quickly and efficiently bring all composite materials to a natural luster.



**Achieve a lustrous,
natural shine in just seconds!**



Pre-Polish



High Shine

By simply using a low speed with light pressure, you can achieve a beautiful, natural shine in a matter of seconds. No paste or water necessary.

The COMPOSITE SINGLES are single-use polishers designed to quickly produce a brilliant luster on all composite materials in the market. This new line of polishers allows for ergonomic access to all surfaces of the tooth for a predictable, high-quality restoration.

Intended for use with all hybrid composites, micro-filled composite restoratives, and high-density restorative materials.

- Beautifully Polishes All Composites Materials
- Single-Patient Use
- Saves Time from Processing
- Simple Two-Step System

See Reverse Side For Ordering Options

MEISINGER USA

10150 E. Easter Avenue • Centennial, CO 80112 • USA

Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407

customerservice@meisingerusa.com • www.meisingerusa.com

   @meisingerusa

Meisinger since 1888

Composite SINGLES

Scan to learn more about the Composite SINGLES on meisingerusa.com



50 Burs Per Pack • RA Shank

Single-patient use, two-step polishing of all composite restorations.

9754M 9754F



50	Grit	Pre-Polish	High Shine
Diameter (1/10mm)		055	055
Length (mm)		8.7	8.7
RA Shank		9754M-055-RA	9754F-055-RA

9755M 9755F



50	Grit	Pre-Polish	High Shine
Diameter (1/10mm)		040	040
Length (mm)		6.9	6.9
RA Shank		9755M-040-RA	9755F-040-RA

9756M 9756F



50	Grit	Pre-Polish	High Shine
Diameter (1/10mm)		060	060
Length (mm)		6.4	6.4
RA Shank		9756M-060-RA	9756F-060-RA

9757M 9757F



50	Grit	Pre-Polish	High Shine
Diameter (1/10mm)		110	110
Length (mm)		1.6	1.6
RA Shank		9757M-110-RA	9757F-110-RA

NEW TO THE COMPOSITE SINGLES?

This assorted pack is perfect for you!



#COSII

Composite SINGLES Assortment Pack

This kit is the perfect way to get started with the Composite SINGLES. It includes an assortment of two of every shape and grit.

A Perfect Pair: Twist Finisher Carbides

Scan to learn more about the Twist Finishers on meisingerusa.com




close-up of the unique
helical fluting

For Superior Composite Trimming and Finishing

- Provides an extremely smooth surface with virtually no surface striations
- Helical flute design ensures continuous contact, eliminating chatter
- Complements all composite materials
- Made in Germany

5 Burs Per Pack FG Shank



5	Figure	D0132	D0133	D0134	D0135	D0379	D0379
Diameter (1/10mm)		008	010	014	014	018	023
Fine (12 Flute)		D0132-008-FG	D0133-010-FG	D0134-014-FG	D0135-014-FG	D0379-018-FG	D0379-023-FG
Extra-Fine (20 Flute)		D0132F-008-FG	D0133F-010-FG	D0134F-014-FG	D0135F-014-FG	D0379F-018-FG	D0379F-023-FG
Ultra-Fine (30 Flute)		D0132U-008-FG	D0133U-010-FG	D0134U-014-FG	D0135U-014-FG	D0379U-018-FG	D0379U-023-FG

MEISINGER USA

10150 E. Easter Avenue • Centennial, CO 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
customerservice@meisingerusa.com • www.meisingerusa.com

@meisingerusa



Contact your local dealer
representative to order today.