

# LUSTER® Paste

## High-gloss on your ceramic restorations

The MEISINGER LUSTER® paste is a 1-step diamond polishing paste and is used for high-gloss polishing of dental ceramic materials. The high-quality LUSTER® paste contains natural diamond particles with a grain size of 8 - 12 µm (fine) and 16 - 18 µm (medium). The paste gives perfect high gloss to all ceramic restorations and is suitable for both extraoral and intraoral polishing. This simplifies the polishing of occlusal surfaces and the maintenance of existing restorations.

Made in  
Germany



### At a glance

- Effective polishing thanks to natural diamond particles
- Perfect high gloss results on all ceramic restorations
- Suitable also for intraoral application
- Economical dosage

Our sales force would be pleased to  
answer your questions: +49 2131 2012-303

Meisinger  
since  
1888

# LUSTER<sup>®</sup> Paste

|   | Figure | Description                          | Grit                  | Filling Amount    | Amount  |
|---|--------|--------------------------------------|-----------------------|-------------------|--|
|  | DP2    | Polishing paste, fine<br>Container   | 8 – 12 $\mu\text{m}$  | 3 g<br>Container  |  |
|  | DP4    | Polishing paste, medium<br>Container | 16 – 18 $\mu\text{m}$ | 3 g<br>Container  |  |
|  | DPS2   | Polishing paste, fine<br>Applicator  | 8 – 12 $\mu\text{m}$  | 2 g<br>Applikator |  |

CE

Available from the dental trade.

|  |   |
|--|---|
| Your trade depot:<br><br>_____<br><br>_____<br><br>_____ | Your delivery address:<br><br>_____<br><br>_____<br><br>_____ |
|--|---|

  
**Hager & Meisinger GmbH**  
 Hansemannstr. 10 • 41468 Neuss • Germany  
 Phone: +49 2131 2012-0  
 info@meisinger.de • www.meisinger.de

**Meisinger USA, L.L.C.**  
 10200 E. Easter Avenue • Centennial, Colorado 80112 • USA  
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407  
 Toll free: +1 (866) 634-7464  
 info@meisingerusa.com • www.meisingerusa.com

85FL123EN - 0817