



One Size Doesn't Fit All...

smaller patients need smaller diamonds!



MEISINGER's Micro-Diamonds are the best fit for your pediatric patients. With particularly small diameters and a precise construction, these Micro-Diamonds are perfect for minimally invasive pediatric procedures.

Micro-Diamonds

- Micro-sizes for maximum preservation of sound tooth structure
- Allow for tighter margins and increased visibility
- Natural diamonds provide efficient cutting for reduced chair time
- Made in Germany

Visit meisingerusa.com or see reverse side for ordering options

MEISINGER USA

10150 E. Easter Avenue • Centennial, CO 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
customerservice@meisingerusa.com • www.meisingerusa.com

   @meisingerusa



5 Burs Per Pack • FG Shank

- Ideal for narrow working areas and tight margins
- Built with premium, natural diamonds for superior cutting
- Micro sizes for maximum conservation
- Made in Germany



Figure #	M801	M833	M838	M880	M889	M890	M895	M896
Diameter (1/10mm)	007	009	007	007	007	008	007	007
Length (mm)	-	1.1	2.1	3.6	2.1	3.6	2.1	3.6
FG	M801-007-FG	M833-009-FG	M838-007-FG	M880-007-FG	M889-007-FG	M890-008-FG	M895-007-FG	M896-007-FG



Micro-Diamond Assorted Kit

This kit is the perfect way to start your collection of our Micro-Diamonds, as well as a great organization tool for your pediatric practice.


Popular in Pediatric Offices:

Scan to learn more about our Carbides on meisingerusa.com




Carbides designed for successful pediatric procedures

- Made with premium, tungsten carbide to provide efficient cutting and superior strength
- Diameters smaller than 1 mm to ensure safe cutting in narrow areas
- Various shapes available for specific procedures
- Made in Germany



US#	1/4	1/2	-	1	-
Figure #	HM1	HM1	HM1	HM1	HM1
Diameter (1/10mm)	005	006	007	008	009
FG	HM1-005-FG	HM1-006-FG	HM1-007-FG	HM1-008-FG	HM1-009-FG
FGXL (SU)	HM1-005-SU	HM1-006-SU	-	HM1-008-SU	-
RA	HM1-005-RA	HM1-006-RA	HM1-007-RA	HM1-008-RA	HM1-009-RA



US#	245
Figure #	HMUN245
Diameter (1/10mm)	009
Length (mm)	2.9
FG	HMUN245-009-FG
SS	HMUN245-009-SS



US#	329	330	330L
Figure #	HM7	HM7	HM7L
Diameter (1/10mm)	006	008	008
Length (mm)	1.7	1.8	3.6
FG	HM7-006-FG	HM7-008-FG	HM7L-008-FG
RA	-	HM7-008-RA	-
SS	-	HM7-008-SS	-

Visit meisingerusa.com to view our entire collection of Carbide Burs

MEISINGER USA

10150 E. Easter Avenue • Centennial, CO 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
customerservice@meisingerusa.com • www.meisingerusa.com



Contact your local MEISINGER Representative to order today.