


# Twist Polishers by MEISINGER

## Efficient high-gloss polishing of composites and ceramics

LUSTER® Twist Polisher Kits are especially developed for high-gloss polishing of all ceramic and composite restorations. The flexible polishing lamellae make it possible to adapt to any surface and render flawless results in areas that are difficult to reach, such as occlusal surfaces.



**LUS06**  
**LUSTER® Intraoral Twist Kit**  
Intraoral high-gloss polishing of all ceramic restorations



**LUS05**  
**LUSTER® Extraoral Twist Kit**  
Extraoral high-gloss polishing of all ceramic restorations



**2663**  
**Twist Polishing Kit**  
Intraoral high-gloss polishing of all composite restorations



### At a glance

- Flexible polishing lamellae can adapt to any surface structure
- Specifically for areas that are difficult to reach, such as occlusal surfaces
- Low heat production due to permanent air cooling of the material
- Designed for gentle action, which preserves the original surface structure
- Graduated diamond grit system allows for short processing times and a long service life
- User friendly and intuitive handling

**MEISINGER USA, L.L.C.**

10200 E. Easter Ave. • Centennial • Colorado • 80112

Tel.: +1 (303) 268-5400 • Fax: +1 (303) 268-5407







www.meisingerusa.com • customerservice@meisingerusa.com

*Meisinger* since 1888

# Twist Polishers by MEISINGER

Fax for ordering (303) 268-5407

	Fig.	Size	Shank	Pieces/ Kit	Reorder Qty		recomm. / max. rpm	Qty 
	LUS05	-	-	1	1	LUSTER® Extraoral Twist Kit	-	
	9771G	170	HP (104)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771M	170	HP (104)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771F	170	HP (104)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771C	170	HP (104)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	LUS06	-	-	1	1	LUSTER® Intraoral Twist Kit	-	
	9771G	140	RA (204)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771M	140	RA (204)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771F	140	RA (204)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	9771C	140	RA (204)	1	2	ceramic restorations	10,000 - 12,000 / 20,000 rpm	
	2663	-	-	1	1	Twist Polishing Kit	-	
	9769M	140	RA (204)	2	2	composite restorations	10,000 - 12,000 / 20,000 rpm	
	9769M	100	RA (204)	2	2	composite restorations	10,000 - 12,000 / 20,000 rpm	
	9769F	140	RA (204)	2	2	composite restorations	10,000 - 12,000 / 20,000 rpm	
	9769F	100	RA (204)	2	2	composite restorations	10,000 - 12,000 / 20,000 rpm	

Diamond polishing system:  Pre-polishing  Polishing  High shine polishing  Mirror finish polishing  
 Pre-polishing  High shine polishing

CE 0044

Available from your local dental dealer.

Your dental dealer:

Your delivery address:

---



---



---



---



---



---



---



---

**Hager & Meisinger GmbH**

Hansemannstr. 10 • 41468 Neuss • Germany  
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222  
 info@meisinger.de • www.meisinger.de

**Meisinger USA, L.L.C.**

10200 E. Easter Avenue • Centennial, Colorado 80112 • USA  
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407  
 Toll free: +1 (866) 634-7464  
 info@meisingerusa.com • www.meisingerusa.com